NanoZoomer[®] S60v2MD Slide scanner system

C16600-21MDEU

Scan various size slides with this regulatory compliant model



Multi-slide Type Scanning

Up to 60 standard-sized slides or 30 mega-sized slides

Selectable Scanning Modes

Fully and semi-Automatic scanning available

Regulatory Compliant

IVDR (EU) UKCA (UK) IVDO (CH)

Highspeed Scanning

Lowoperational Workload

60 s (20× mode) 75 s (40× mode) * When it scans an area of 15 mm × 15 mm square with 5 focus points.


Assistant for image quality check



Flexible multi-slide scanning options at constant high-speed

Multi-slide type handling

The NanoZoomer® S60v2MD has three slots available, which allows to scan both standard and mega-sized slides in one batch.



Users can use up to three slots using standard or mega-sized slides interchangeably based on their preferences.



More productive and convenient

- Profile functions
 - Scan profile

Save scan conditions such as scan area and sample size threshold for tissue recognition.



• Output profile You can choose multiple files for saving an image.



High-speed scanning

Achieve these outstanding scan timings on both standard and mega-sized slides.



Improved scanning workflow solutions

Choose your preferred scanning mode.

Fully-Automatic scanning



All scanning processes work automatically.

Maintain an optimized system condition



Optimum image quality and color balance are maintained by

the operational software "NZAcquireMD".


Quality check

QC (Quality Check) mode is available to allow users to check image quality before finalizing Whole Slide Imaging.





These scan timings can be achieved when it scans an area of 15 mm × 15 mm square with 5 focus points.



Semi-Automatic scanning



Option to set-up scanning conditions such as the scan area or resolution and to assign profiles for each slide.



.....




Checking

System configuration

Software Software Hardware NZAcquireMD NZViewMD Scanner (For operation) (For viewing)

Dimensional outlines (Unit: mm)



Specifications

•			
Product name			NanoZoomer [®] S60v2MD Slide scanner system
Product number			C16600-21MDEU
Scanning	20× mode		Approx. 60 s
speed *1	40× mode		Approx. 75 s
Objective lens			20× N.A. 0.75
			User can select 20× or 40× mode
			at start of scanning
Compatible glass slides			25.0 mm to 26.0 mm × 75.0 mm to 76.0 mm (Thickness 0.9 mm to 1.2 mm)
			51.0 mm to 52.0 mm × 75.0 mm to 76.0 mm
			(Thickness 0.9 mm to 1.2 mm)
Slide loader	Standard-sized slide *2		Up to 60 slides (20 slides × 3 cassettes) *3
	Mega-sized slide *2		Up to 30 slides (10 slides × 3 cassettes) *3
Scanning resolution	20× mode		Approx. 0.46 µm/pixel
	40× mode		Approx. 0.23 µm/pixel
Focusing method			Pre-Focus map
Z-stack feature			Included
Image compression			JPEG compression
		1D Barcodes	Code 39, Code 128, Interleaved 2 of 5, Codabar, EAN-8 and UPC-E
Readable barcodes			DataMatrix (ECC200)
		2D Barcodes	QR code (QR Code Model-1) QR Code Model-2)
Power supply			AC 100 V to AC 240 V
Power consumption (Scanner only)			Approx. 120 VA

*1 When it scans an area of 15 mm × 15 mm square with 5 focus points

*2 Standard-sized is a single slide, mega-sized is a double-width, large slide. *3 Up to 3 cassettes can be set with combination of standard and mega-sized slides.

Intended Use

NanoZoomer® S60v2MD Slide scanner system ("NanoZoomer® System") is an automated digital slide creation, viewing, and management system. The NanoZoomer® System is intended for in vitro diagnostic use as an aid to the pathologist to review and interpret digital images of surgical pathology slides prepared from formalin-fixed paraffin-embedded ("FFPE") tissue. The NanoZoomer® System is not intended for use with frozen section, cytology, or non-FFPE hematopathology specimens

The NanoZoomer® System comprises the Scanner and the NZViewMD Software. The NanoZoomer® System is for creation and viewing of digital images of scanned glass slides that would otherwise be appropriate for manual visualization by conventional light microscopy. It is the responsibility of a qualified pathologist to employ appropriate procedures and safeguards to assure the validity of the interpretation of images obtained using NanoZoomer® System.

• NanoZoomer is a registered trademark of Hamamatsu Photonics K.K. (EU, Japan, UK, USA)

- The product and software package names noted in this brochure are trademarks or registered trademarks of their respective manufacturers.
- Subject to local technical requirements and regulations, availability of products included in this brochure may vary. Please consult with your local sales representative. • The product described in this brochure is designed to meet the written specifications, when used strictly in accordance with all instructions.

© 2022 Hamamatsu Photonics K.K.

HAMAMATSU PHOTONICS K.K. www.hamamatsu.com

Manufacturer

HAMAMATSU PHOTONICS K.K., Systems Division Joko Factory

812 Joko-cho, Higashi-ku, Hamamatsu-City, Shizuoka-Pref. 431-3196, Japan Telephone: (81)53-431-0124, Fax: (81)53-433-8031

E-mail: export@sys.hpk.co.jp

Authorised representative

HAMAMATSU PHOTONICS DEUTSCHLAND GMBH

Arzbergerstr. 10, 82211 Herrsching am Ammersee, Germany EC REP E-mail: pms-med@hamamatsu.eu

HAMAMATSU PHOTONICS UK LIMITED

UK REP

2 Howard Court, 10 Tewin Road, Welwyn Garden City, Hertfordshire, AL7 1BW, UK

E-mail: pms-med@hamamatsu.co.uk



Importers

HAMAMATSU PHOTONICS DEUTSCHLAND GMBH

Arzbergerstr. 10, 82211 Herrsching am Ammersee, Germany Telephone: (49)8152-375-0, Fax: (49)8152-265-8 E-mail: info@hamamatsu.de

HAMAMATSU PHOTONICS FRANCE S.A.R.L.

19 Rue du Saule Trapu, Parc du Moulin de Massy, 91882 Massy Cedex, France Telephone: (33)1 69 53 71 00, Fax: (33)1 69 53 71 10 E-mail: infos@hamamatsu.fr

HAMAMATSU PHOTONICS NORDEN AB

Torshamnsgatan 35 16440 Kista, Sweden Telephone: (46)8-509 031 00, Fax: (46)8-509 031 01 E-mail: info@hamamatsu.se

HAMAMATSU PHOTONICS ITALIA S.R.L.

Strada della Moia, 1 int. 6, 20044 Arese (Milano), Italy Telephone: (39)02-93 58 17 33, Fax: (39)02-93 58 17 41 E-mail: info@hamamatsu.it