CAMERA LINE UP CATALOG




















We have a diverse lineup of cameras that support a wide range of wavelengths from X-rays to the near-infrared and support a variety of applications.

Life science

Super resolution microscopy


Cells are observed with higher spatial resolution than diffraction limit by the super resolution microscopy.



Camera: ORCA®-Quest Super resolution imaging system: VT-iSIM
Data courtery of: Steven Coleman (Visitech international Ltd.)

Bioluminescence measurements


Ligand-stimulated binding of NanoLuc-Arrestin to GPCR-mVenus and its intracellular uptake are observed by simultaneous two-wavelength luminescence imaging.



Data courtery of: Masataka Yanagawa (Tohoku university)

Light sheet microscope


The zebrafish larvae brain function during its natural behavior is observed with light sheet fluorescence microscope.



Data courtery of: Drew Robson (Max Planck Institute for Biological Cybernetics)

Observation of cultured cells

Cells cultured in one well of a microplate are observed by high-resolution imaging with fluorescent images




* Displayed with pseudo color by image processing.

Quantum Technology

Quantum computing (Neutral atom, Ion)

Position and quantum state of Rb atoms, trapped and arranged one by one in vacuum, are diagnosed via fluorescence.



Camera: ORCA®-Quest Data courtery of: Takashi Yamamoto, Toshiki Kobayashi (Osaka university)

Quantum optics

qCMOS® camera is used for absorption imaging with quantum light source to compare between 1 photon event and 2 photon event images.

1 photon event

2 photon event

Data courtery of: Miles Padgett (University of Glasgow)

Astronomy

Lucky imaging

Wide field of view and low-noise imaging is used to obtain a clear image of the stars by integrating, from among many acquired images, that are less affected by atmospheric turbulence.



* Displayed with pseudo color by image processing.


Foreign object detection

Food inspection

Synchrotron imaging

X-ray phase contrast CT image of mouse embryo

The mouse embryo is observed using the synchrotron X-ray.




ORCA®_Quest combined with High resolution X-ray imaging system (M11427) Data courtery of: SPring-8 BL20B2 beamline by Dr. Masato Hoshino, Senior researcher in Japan Synchrotron Radiation Research Institute (JASRI)

Semiconductor inspection


Transmission observation of Si wafer

The pattern formed on the backside of the Si wafer is observed transmitted through the front side by infrared imaging.



Semiconductor device observation

The pattern under the Si layer is observed by infrared imaging.




Structure observation of semiconductor devices

The interior structure of a semiconductor device is analyzed at the nano-level by high-resolution imaging using an electron microscope.



visible light is detected by the infrared imaging.


Small stones mixed in coffee beans that are difficult to see with



Analysis / Spectroscopy

Material identification

Infrared imaging identifies materials that are difficult to distinguish in visible light, such as PVC, acrylic, PET, and PS.





Imaged at multiple wavelengths


* Displayed with pseudo color by image processing.

Wavelength range	Visible to near-infrared						Near-infrared						
Name	ORCA®-Quest 2 qCMOS® camera	ORCA®-Fire Digital CMOS camera	ORCA®-Fusion BT Digital CMOS camera	ORCA®-Fusion Digital CMOS camera		ORCA®-Flash4.0 LT3 Digital CMOS camera	ORCA®-spark Digital CMOS camera	TDI camera		InGaAs	s camera		InGaAs line scan camera
Туре	C15550-22UP	C16240-20UP	C15440-20UP	C14440-20UP	C13440-20CU	C11440-42U40	C11440-36U	C10000-801	C16741-40U	C14041-10U	C12741-03	C12741-11	C15333-10E04
Appearance		0						0		(0)			0
Image sensor type				Area sensor				TDI sensor		Area	sensor		Line sensor
Sensitivity wavelength range (nm)	250 to 1000	250 to 1000	350 t			1000		200 to 1000	400 to 1700	950 t	o 1700	900 to 1550	950 to 1700
(Spectral response: See P5)			3	—• —		6 ——	 6		*3		<u>8</u> ——	—9 —	
Effective number of pixels (H × V)	4096 × 2304	4432 × 2368	2304 × 2304	2304 × 2304	2048 × 2048	2048 × 2048	1920 × 1200	2048 × 128	1280 × 1024	320 × 256	640 × 512	640 × 512	1024 × 1
Pixel size ((H) μm × (V) μm)	4.6 × 4.6	4.6 × 4.6	6.5 × 6.5	6.5 × 6.5	6.5 × 6.5	6.5 × 6.5	5.86 × 5.86	12 × 12	5 × 5		20 × 20		12.5 × 12.5
Effective area ((H) mm × (V) mm)	18.841 × 10.598	20.387 × 10.892	14.976 × 14.976	14.976 × 14.976	13.312 × 13.312	13.312 × 13.312	11.25 × 7.03	24.58 × 1.536	6.40 × 5.12	6.4 × 5.12	12.8 × 10.24	12.8 × 10.24	12.8 × 0.0125
Full well capacity (electrons) typ.*1	7000	20 000	15 000	15 000	30 000	30 000	33 000	80 000	-		-	300 000	-
Dynamic range typ.*1	23 000:1	20 000:1	21 400:1	21 400:1	37 000:1	33 000:1	5000:1	1600:1	-		-	-	-
Cooling method	Forced-air cooled/ Water cooled	Forced-air cooled	Forced-air cooled/ Water cooled	Forced-air cooled/ Water cooled	Forced-air cooled/ Water cooled	Forces-air cooled	-	-	Forced-air cooled/ Natural-air cooled	Forced-	air cooled	Forced-air cooled/ Water cooled	-
Cooling temperature (°C)*1	-35 (Water cooled)	+20	-15 (Water cooled)	-15 (Water cooled)	-30 (Water cooled)	+10	-	-	+15 (Forced-air cooled)	4	·10	-70 (Water cooled)	-
Readout speed (frame/s) (Full resolution)*1	120	115	89.1	89.1	100	30	64.9	50 kHz (Line rate)	71.53	216.6	59.774	7.2	40 kHz (Line rate)
Readout noise (electrons) rms typ.*1	0.30	1.0	0.7	0.7	1.4	1.5	6.6	50	-		-	500	-
Dark current (electrons/pixel/s) typ.*1	0.006 (Water cooled)	0.6	0.7 (Water cooled)	0.2 (Water cooled)	0.006 (Water cooled)	0.6	-	-	-		-	130 (Water cooled)	-
Interface	CoaXPress (Quad CXP-6)/ USB 3.1 Gen1	CoaXPress (Quad CXP-6)/ USB 3.1 Gen1	CoaXPress (Dual CXP-6)/ USB 3.0 *2	CoaXPress (Dual CXP-6)/ USB 3.0 *2	Camera Link/USB 3.0 *2	USB 3.1 Gen 1	USB 3.0 *2	Camera Link	USB 3.1 Gen 1	USB 3.0 *2	USB 3.0 *2/EIA	Camera Link	Gigabit Ethernet
Applications	Life science imaging Quantum technology Astronomy Semiconductor inspection Synchrotron imaging Electron microscope	Life science imaging Synchrotron imaging Electron microscope	Life science imaging Synchrotron imaging Electron microscope	Life science imaging Semiconductor inspection Synchrotron imaging Electron microscope	Life science imaging Semiconductor inspection Synchrotron imaging Electron microscope	Life science imaging Semiconductor inspection	Life science imaging Synchrotron imaging	Life science imaging Semiconductor inspection	Life science imaging Semiconductor inspection	Semiconductor inspection Food inspection Analysis/spectroscopy		Life science imaging Semiconductor inspection	Semiconductor inspection

Camera type			Board type camera for OEM					
Name	Scientific CMOS I	ooard level camera	Digit	TDI board level camera				
Туре	C11440-62U C11440-52U30		C13949-50U C13770-50U		C13752-50U	C10000-A01		
Appearance	0		0	0	0	9		
Image sensor type	Area sensor			TDI sensor				
Sensitivity wavelength range (nm)	350 to	1000		200 to 1000				
(Spectral response: See P5)	$\overline{}$							
Effective number of pixels (H × V)	2048 × 2048		4096 × 3008	2464 × 2056	2048 × 1544	2048 × 128		
Pixel size ((H) μm × (V) μm)	6.5 × 6.5			12 × 12				
Effective area ((H) mm × (V) mm)	13.312 × 13.312		14.13 × 10.37	8.50 × 7.09	7.06 × 5.32	24.53 × 1.536		
Full well capacity (electrons) typ.*1	30 000			80 000				
Dynamic range typ.*1	20 000:1 18 000:1			1600:1				
Readout speed (frame/s) (Full resolution)*1	30		15	40	65	50 kHz (Line rate)		
Readout noise (electrons) rms typ.*1	2.1	2.3	2.3		50			
Interface	USB	3.0 *2		Camera Link				
Applications	Contact us	Contact us	Contact us	Contact us	Contact us	Contact us		


^{*1} Depends on the mode and conditions. For details, refer to each product catalog. *2 Equivalation to USB 3.1 Gen1

For X-ray				
X-ray sCMOS camera	ORCA®-Lightning X X-ray sCMOS camera			
C12849-111U	C15606-101P C15606-102			
0	0 1			
Area sensor	Area sensor			
25 kV to 90 kV (Recommended X-ray tube voltage range)	25 kV to 70 kV (Recommended X-ray tube voltage range)			
2048 × 2048	4608 × 2592			
6.5 × 6.5	5.5 × 5.5			
13.312 × 13.312	25.344 × 14.256			
30 000	38 000			
18 000:1	15 000:1			
30	30			
2.3	3.0			
USB 3.0 *2	CoaXPress (Quad CXP-6)			
Synchrotron imaging	Synchrotron imaging			



OEM CAMERA


We design and manufacture OEM cameras specific to each customer. We provide various types of cameras with options such as shape, sensor, interface, cooling method, software, etc. to meet customers' requests. The measurement wavelength range covers not only the visible range but spans widely from X-ray to infrared.



Shorten delivery time with simulation technology

We can perform imaging simulations that match the characteristics of various cameras (wavelength, sensitivity, speed, etc.). By using this technology, we can shorten the process of repeating design and trial production, and provide cameras that meet your purpose efficiently and in a short time.

Flow from inquiry to delivery




You can try the simulation on our website. Access it from the URL below.

Camera Simulation Engine URL: https://camera.hamamatsu.com/all/en/learn/camera_simulation_engine.html

SOFTWARE

We provide a common camera library "DCAM-API®," Hamamatsu Photonics software that can maximize the characteristics of your camera, and a tool "DCAM-SDK," that allows you to build your own control software. Through DCAM-API®, even if the camera or interface is changed, the software modification/change can be minimized.



Third-party software

Our cameras can be controlled by combining our cameras and peripherals with software from each microscope manufacturer, bioimaging software, or the following software.

Plugins that are compatible with third-party software

Software	Manufacturer	os
LabView	National Instruments	Windows®
MATLAB	The MathWorks	Windows®
μManager	Open source	Windows [®]
EPICS	Open source	Linux®

*For details on external software, please contact the manufacturer

Please download plugins from the URL below. URL: https://dcam-api.com/third-party-plugins/

For details, please refer to the following link.

https://www.hamamatsu.com/all/en/product/cameras/driver-software.html

RELATED PRODUCTS

Imaging optical system

We also have a lineup of Imaging optical systems to expand usability of our cameras, such as multi wavelength imaging and High resolution X-ray imaging system.



Image splitting optics W-VIEW GEMINI A12801

Product details page URL: https://www.hamamatsu.com/all/en/product/optical-components/imagesplitting-optics.html



High resolution X-ray imaging system M11427

Product details page URL: https://www.hamamatsu.com/all/en/product/cameras/high-resolution-x-rayimaging-system.html

X-ray line scan camera/X-ray TDI camera


We have a lineup of X-ray non-destructive inspection cameras that can be used in-line. Since it is possible to inspect the inside of substances that cannot be seen with visible light or infrared light, these cameras are suitable for foreign matter inspection of foods and pharmaceuticals, defect inspection of printed circuit board, etc.



X-ray line scan camera C14300 series



X-ray TDI camera C12300 series



For details, please refer to the following link.

https://www.hamamatsu.com/all/en/applications/non-destructive-testing.html

- DCAM-API is registered trademark of Hamamatsu Photonics K.K. (EU, Japan, UK, USA).
- ORCA and qCMOS are registered trademarks of Hamamatsu Photonics K.K. (China, EU, Japan, UK, USA).
- Windows is a registered trademark of Microsoft Corporation in the USA (and other countries)
- The product and software package names noted in this brochure are trademarks or registered trademarks of their respective manufacturers.
- Subject to local technical requirements and regulations, availability of products included in this brochure may vary. Please consult your local sales representative. • The products described in this brochure are designed to meet the written specifications, when used strictly in accordance with all instructions
- The spectral response specified in this brochure is typical value and not guaranteed.
- Specifications and external appearance are subject to change without notice.
- © 2024 Hamamatsu Photonics K.K.

HAMAMATSU PHOTONICS K.K. www.hamamatsu.com

812 Joko-cho, Chuo-ku, Hamamatsu City, 431-3196, Japan, Telephone: (81)53-431-0124, Fax: (81)53-433-8031, E-mail: export@sys.hpk.co.jp

S12 JOKO-CIO, CHUO-KU, Harnamatsu City, 431-3196, Japan, Telephone: (31) 33-431-U124, Fax: (81) 35-433-8031, E-mail: export@sys.npk.co.jp

U.S.A.: HAMAMATSU CORPORATION: 360 Foothill Road, Bridgewater, NJ 08807, U.S.A., Telephone: (1)908-231-0960, Fax: (1)908-231-1218

Germany: HAMAMATSU PHOTONICS DEUTSCHLAND GMBH: Arzbergerstr. 10, 82211 Herrsching am Ammersee, Germany, Telephone: (49)8152-375-0, Fax: (49)8152-265-8 E-mail: info@hamamatsu.de

France: HAMAMATSU PHOTONICS FRANCE S.A.R.L.: 19 Rue du Saule Trapu, Parc du Moulin de Massy, 91882 Massy Cedex, France, Telephone: (33) 1 69 53 71 00. Fax: (34)8-508 1 70. Telephone: (34)8-508 1 70. Fax: (34)8-508 1 70. Fax: (35)8-508 1 70. Fax: (35)8-508 1 70. Fax: (36)8-508 1 70. F