

Fast kinetics enzyme assay High-throughput screening of prolyl isomerase inhibitors

Abstract

The peptidyl-prolyl cis-trans isomerases (PPIases), or rotamases, are a family of ubiquitous enzymes that catalyze the cis-trans isomerization of the peptide bond on the N-terminal side of proline residues in proteins. PPIases are divided into three classes, cyclophilins (CyPs), FK-506 binding proteins (FKBPs) and the Pin1/parvulin class.

Cyclophilin(Cyp) is one of the key enzymes and plays critical roles during many biological processes and the unique enzyme Pin1 plays a key role in the pathogenesis of Alzheimer's disease and many human cancers. Pin1 also regulates the cell cycle and is a necessary enzyme for cell division. For these reasons, cyclophilins have emerged as potential drug targets for several diseases. Therefore, it is extremely important to screen for novel small molecule cyclophilin inhibitors.

Recently, a real time fluorescence monitoring system for Cyclophilin A was reported on a high-throughput screening format. PPIase activity has been assessed by assays based on protease-coupled isomer specific proteolysis using tetrapeptide derivatives, Ala-Xaa-Pro-Yaa-(pNA/MCA). Cyclosporin A inhibits the PPIase activity, and it has been used as tool compound during the setup.



Filter set is available.

Please contact Hamamatsu representative with the part number A10343-31 for AMC filter set.

Suc-Ala-Ala-Pro(cis)-Phe-MCA

↓ CypA

Suc-Ala-Ala-Pro(trans)-Phe-MCA

↓ chymotrypsin

Suc-Ala-Ala-Pro-Phe-COOH + AMC

Screening using FDSS7000/6000/3000/uCell

AMC (λ_{ex} =340 nm to 380 nm; λ_{em} =460 nm)




Fig. 1: Cyclophilin A (CypA) peptidyl-prolyl isomerase (PPIase) assay.

The protease-couple CypA PPIase assay carried out by measuring the increase of fluorescence intensity from liberated 7-amino-4-methylcoumarin(AMC). Then 15 μ M fluorogenic peptide substrate that had been proline isomerized was hydrolyzed by 50 μ M chymotrypsin, an isomer-specific protease.

(A) Time course of CypA activity. The reaction was performed for 60s, and at the beginning of 10 sec, the reaction was shown because the CypA-catalyzed cis-trans isomerization rate is extremely fast.

The fluorescence of the liberated AMC in each well was detected simultaneously by the FDSS6000 real time fluorescence detector. Activity was measured at the indicated concentrations of enzyme.

(C) Concentration dependence of CsA inhibition of CypA activity. Concentration of CypA and CsA are expressed as follows.

From T. Mori, et al., J. Biomol. Screen. 14(4), 419-424, copyright© 2009 by SAGE Publications. Reprinted by permission of SAGE Publications.




Fig. 2:

Evaluation of cyclophilin A (CypA) peptidyl-prolyl cis-trans isomerase (PPIases) assay performance. Forty 384-well plates that contain 16 control (CypA 7nM) and blank (CypA none) wells, respectively, were used to determine the suitability of the assay for high-throughput screening. (A) The conversion ratio of substrate for each well was recorded, and average values corresponding to each plate were plotted. (B) The Z' factors were calculated for the each plate.

From T. Mori, et al., J. Biomol. Screen. 14(4), 419-424, copyright© 2009 by SAGE Publications. Reprinted by permission of SAGE Publications.

Summary

Dr. Mori and Prof. Uchida (Tohoku university) have established the High-Throughput Screening method for Prolyl Isomerase Inhibitors using a Real-Time Fluorescence Monitoring System. To use this method, Their group have finished an HTS successfully.

References

1. Mori, T. et al. A dual inhibitor against prolyl isomerase Pin1 and cyclophilin discovered by a novel real-time fluorescence detection method. *Biochem. Biophys. Res. Commun.* 406, 439-443 (2011).
2. Mori, T. et al. Use of a real-time fluorescence monitoring system for high-throughput screening for prolyl isomerase inhibitors. *J. Biomol. Screen.* 14, 419-424 (2009).
3. Tataru, Y. et al. Dipentamethylene thiuram monosulfide is a novel inhibitor of Pin1. *Biochem. Biophys. Res. Commun.* 384, 394-398 (2009).
4. Uchida, T. et al. Identification and characterization of a 14 kDa human protein as a novel parvulin-like peptidyl prolyl cis/trans isomerase. *FEBS Lett.* 446, 278-282 (1999).

★ FDSS is registered trademark of Hamamatsu Photonics K.K. (China, France, Germany, Italy, Japan, U.K., U.S.A.)

★ Product and software package names noted in this documentation are trademarks or registered trademarks of their respective manufacturers.

- Subject to local technical requirements and regulations, availability of products included in this promotional material may vary. Please consult your local sales representative.
- Information furnished by HAMAMATSU is believed to be reliable. However, no responsibility is assumed for possible inaccuracies or omissions. Specifications and external appearance are subject to change without notice.

© 2013 Hamamatsu Photonics K.K.

HAMAMATSU PHOTONICS K.K. www.hamamatsu.com

HAMAMATSU PHOTONICS K.K., Systems Division

812 Joko-cho, Higashi-ku, Hamamatsu City, 431-3196, Japan, Telephone: (81)53-431-0124, Fax: (81)53-435-1574, E-mail: export@sys.hpk.co.jp

U.S.A.: Hamamatsu Corporation: 360 Foothill Road, Bridgewater, N.J. 08807-0910, U.S.A., Telephone: (1)908-231-0960, Fax: (1)908-231-1218 E-mail: usa@hamamatsu.com

Germany: Hamamatsu Photonics Deutschland GmbH: Arzbergerstr. 10, D-82211 Herrsching am Ammersee, Germany, Telephone: (49)8152-375-0, Fax: (49)8152-2658 E-mail: info@hamamatsu.de

France: Hamamatsu Photonics France S.A.R.L.: 19, Rue du Saule Trapu, Parc du Moulin de Massy, 91882 Massy Cedex, France, Telephone: (33)1 69 53 71 00, Fax: (33)1 69 53 71 10 E-mail: infos@hamamatsu.fr

United Kingdom: Hamamatsu Photonics UK Limited: 2 Howard Court, 10 Tewin Road, Welwyn Garden City, Hertfordshire AL7 1BW, United Kingdom, Telephone: (44)1707-294888, Fax: (44)1707-325777 E-mail: info@hamamatsu.co.uk

North Europe: Hamamatsu Photonics Norden AB: Torshamnsgatan 35 SE-164 40 Kista, Sweden, Telephone: (46)8-509-031-00, Fax: (46)8-509-031-01 E-mail: info@hamamatsu.se

Italy: Hamamatsu Photonics Italia S.r.l.: Strada della Moia, 1 int. 6, 20020 Arese (Milano), Italy, Telephone: (39)02-93581733, Fax: (39)02-93581741 E-mail: info@hamamatsu.it

China: Hamamatsu Photonics (China) Co., Ltd.: B1201 JIaming Center, No.27 Dongsanhuan Beilu, Chaoyang District, Beijing 100020, China, Telephone: (86)10-6586-6006, Fax: (86)10-6586-2866 E-mail: hpc@hamamatsu.com.cn

Cat. No. SBIS0085E02
NOV/2013 IP
Created in Japan (PDF)