



# High Speed Data Acquisition Option for FDSS series

Fast acquisition with stable sampling interval allows us to do more. Hamamatsu presents the new "High Speed Data Acquisition Option" for FDSS series; It enables to acquire data at very short sampling time intervals (approx. 10 ms, 96 well plate, 4x4 binning) to capture the fast  $Ca^{2+}$  oscillation and membrane potential of cardiomyocytes, which is suitable for pre-toxicity study in vitro. Not only for cardiomyocyte but for conventional Aequorin  $Ca^{2+}$  assay, high speed gives us different results.




## ES/iPS Cell-derived cardiomyocytes $Ca^{2+}$ oscillation and membrane potential assay



Dose response in beating of human iPSC-derived cardiomyocytes  $Ca^{2+}$ . Measured few minutes after d-Sotalol addition

When measuring calcium ion or membrane potential oscillation in a cardiomyocyte, such like ES or iPSC-derived derived ones, the acquisition speed give us more relevant information. Below is the comparison between when the interval is 115 ms, which is a conventional speed of acquiring data in fluorescence plate imager, and when the interval is 9 ms, which our new function can achieve.




Upper waveform is measured with 115 ms interval, lower waveform is measured in 9 ms interval using the new high speed data acquisition option. Zoomed into the first 10 s measurement.

Waveform becomes different when measured in very short interval time.

FDSS 7000EX

FDSS 77 CELL




Same as the former page, upper waveform is measured with 115 ms interval, lower waveform is measured in 9 msec interval using the new high speed data acquisition option. Zoomed into the first 40 seconds measurement.

Waveform becomes more uniform when measured in very short interval time.

## Aequorin Calcium ion Luminescence Assay

We have tested our high speed data acquisition option in flash luminescence such as aequorin, and found that there is a difference in the peak timing. Measuring flash luminescence with high speed could give us more information.



On the left is the graph measured with 1 s interval, the right is the graph measured with 16 ms interval. Fast acquisition allow us to detect the flash luminescence peak accurately.

**\*\*Note:** Exposure time on the left (1 s interval) were intentionally reduced for the experimental comparison between above two.

★ **FDSS is registered trademark of Hamamatsu Photonics K.K. (China, France, Germany, Italy, Japan, U.K., U.S.A.)**

★ **Product and software package names noted in this documentation are trademarks or registered trademarks of their respective manufacturers.**

- Subject to local technical requirements and regulations, availability of products included in this promotional material may vary. Please consult your local sales representative.
- Information furnished by HAMAMATSU is believed to be reliable. However, no responsibility is assumed for possible inaccuracies or omissions. Specifications and external appearance are subject to change without notice.

© 2013 Hamamatsu Photonics K.K.

**HAMAMATSU PHOTONICS K.K.** [www.hamamatsu.com](http://www.hamamatsu.com)

**HAMAMATSU PHOTONICS K.K., Systems Division**

**812 Joko-cho, Higashi-ku, Hamamatsu City, 431-3196, Japan, Telephone: (81)53-431-0124, Fax: (81)53-435-1574, E-mail: export@sys.hpk.co.jp**

**U.S.A.:** Hamamatsu Corporation, 360 Foothill Road, Bridgewater, N.J. 08807-0910, U.S.A., Telephone: (1)908-231-0960, Fax: (1)908-231-1218 E-mail: usa@hamamatsu.com

**Germany:** Hamamatsu Photonics Deutschland GmbH, Arzbergerstr. 10, D-82211 Hirsching am Ammersee, Germany, Telephone: (49)8152-375-0, Fax: (49)8152-2658 E-mail: info@hamamatsu.de

**France:** Hamamatsu Photonics France S.A.R.L., 19, Rue du Saule Trapu, Parc du Moulin de Massy, 91882 Massy Cedex, France, Telephone: (33)1 69 53 71 00, Fax: (33)1 69 53 71 10 E-mail: info@hamamatsu.fr

**United Kingdom:** Hamamatsu Photonics UK Limited, 2 Howard Court, Welwyn Garden City, Hertfordshire AL7 1BW, United Kingdom, Telephone: (44)1707-294888, Fax: (44)1707-325777 E-mail: info@hamamatsu.co.uk

**North Europe:** Hamamatsu Photonics Norden AB, Torshamnsgatan 35 SE-164 40 Kista, Sweden, Telephone: (46)8-509-031-00, Fax: (46)8-509-031-01 E-mail: info@hamamatsu.se

**Italy:** Hamamatsu Photonics Italia S.r.l., Strada della Moia, 1 int. 6, 20020 Arese (Milano), Italy, Telephone: (39)02-93581733, Fax: (39)02-93581741 E-mail: info@hamamatsu.it

**China:** Hamamatsu Photonics (China) Co., Ltd., B1201 Jiaming Center, No.27 Dongsanhuan Beltu, Chaoyang District, Beijing 100020, China, Telephone: (86)10-6586-6006, Fax: (86)10-6586-2866 E-mail: hpc@hamamatsu.com.cn

Cat. No. SBIS0095E01  
NOV/2013 IP  
Created in Japan (PDF)