



# Spectral Urine Analysis for Non-Invasive, Rapid Diagnosis: The Future of Liquid Biopsy

Early diagnosis plays a pivotal role in the effective treatment of diseases, particularly in oncology. However, the absence of apparent clinical symptoms and the limitations of conventional screening methods present significant challenges. To address this gap, medical science is exploring innovative diagnostic platforms, among which liquid biopsy has emerged as a promising technique. Introduced initially to detect circulating tumor cells (CTCs) in the bloodstream, liquid biopsy now encompasses a broader range of body fluids, including cerebrospinal fluid, saliva, pleural effusion, and most notably, urine.

## Why Urine?

Urine offers a unique window into human health. As a biological fluid that is readily obtainable through non-invasive methods and widely accepted by patients, urine is increasingly being recognized for its diagnostic potential. It reflects the metabolic state of various organs and can indicate both inflammatory and neoplastic processes. With the presence of metabolites and biomarkers, urine becomes an ideal candidate for spectral analysis.

## Traditional Urine Analysis Methods

A cornerstone of diagnostic medicine, urine analysis can be traced back to ancient times. These traditional methods, though still widely used, have limitations in scope, precision, and speed.

### Culture Tests

Used for diagnosing infections, urine cultures are the gold standard for identifying bacterial species and determining antibiotic sensitivity. However they delay diagnosis with incubation times of 24–72 hours.

### Urinalysis Dipstick Tests

Colorimetric tests are commonly used in clinical settings for quick screening. Dipsticks detect a range of parameters such as pH, glucose, protein, ketones, bilirubin, and blood. They are useful, but semi-quantitative, can be subject to user error, and often provide only a basic overview.

### Urine Microscopy

Urine sediment examination under a microscope is used to detect cells (RBCs, WBCs), casts, crystals, and bacteria. Though informative, microscopy is labor-intensive, requires trained personnel, and may lack sensitivity for low-abundance biomarkers.

### Chemical Assays and Enzymatic Methods

These tests measure specific substances (e.g., creatinine, urea, or albumin) with higher accuracy than dipsticks. They aren't always ideal as they require sample preparation, reagents, and laboratory infrastructure, and are typically time-consuming.


# Spectral Techniques in Urine Analysis

Recent technological advancements have enabled the use of optical (spectral) methods to analyze urine at the molecular level, offering rapid, accurate, and non-invasive diagnostic alternatives. These techniques rely on the interaction of light with urine molecules, capturing spectral signatures that reveal both qualitative and quantitative information.

## Raman Spectroscopy

Raman spectroscopy measures the inelastic scattering of monochromatic light—typically from a laser—as it interacts with molecular vibrations. It is highly sensitive and can detect trace amounts of substances. This method has proven especially effective in identifying:

- **Drugs of abuse** (e.g., opioids, amphetamines, MDMA)
- **Therapeutic agents** (e.g., benzodiazepines, antipsychotics)
- **Doping substances** (e.g., erythropoietin, ephedrine)



Spectral characteristics of urine specimens from healthy human volunteers analyzed using Raman chemometric urinalysis (Rametrix), Senger RS, Kavuru V, Sullivan M, Gouldin A, Lundgren S, Merrifield K, et al. (2019) <sup>[1]</sup>



Raman spectra from 235 urine specimens obtained using 735 nm (30 mW) laser excitation. (A) Overlaid vector-normalized spectra. (B) Average urine spectrum (black) with 1 (red) and 2 (blue) standard deviation ranges. Notable peaks appear at 1006  $\text{cm}^{-1}$  (Urea) and 1640  $\text{cm}^{-1}$  (Water).

---

## Fourier Transform Infrared Spectroscopy (FTIR)

FTIR captures high-resolution spectral fingerprints of urine, enabling the detection of proteins, lipids, cells, and even pathogens. Combined with computational techniques such as principal component analysis (PCA) and Soft Independent Modeling of Class Analogy (SIMCA), FTIR has demonstrated strong clinical potential, including the ability to distinguish kidney donors at risk of organ rejection.



Hamamatsu's C16511-01 FTIR engine, high-speed and compact near infrared spectrometer that can be incorporated into inline process.

---

## Visible (VIS) and Near-Infrared (NIR) Spectrometry

These techniques use a broadband light source that transmits or reflects through the urine sample. A spectrometer then detects the specific wavelengths absorbed by the sample. This allows for the assessment of a broad range of biomarkers which can be indicators of several dysfunctions:

- **Hydration and renal function:** Urine osmolality, specific gravity, creatinine
- **Metabolic health:** Uric acid, urea, citrate, oxalate
- **Infection and inflammation:** White/red blood cells, bacteria, albumin, total protein

Although inorganic ions do not directly absorb NIR light, they affect spectral patterns through ion–water interactions. Likewise, urine saturation and crystallization (e.g., calcium oxalate) can be detected optically. Measuring urine saturation and crystallization is vital for preventing, diagnosing, and managing kidney stones and urinary tract conditions. For example, early detection of supersaturated urine can enable timely dietary, hydration, or pharmacological interventions to prevent stone formation.



Hamamatsu's TF series C14486GA minispectrometer, for near-infrared range, compact / thin case with InGaAs linear image sensor mounted.

---

## Fluorescence Spectrometry

This method measures light emissions from substances that have absorbed energy, typically UV or visible light. Fluorescence spectrometry detects endogenous fluorophores such as:

- **NADH and FAD** (indicators of mitochondrial function)
- **Tryptophan** (linked to protein metabolism and cellular activity)
- **Porphyrins** (associated with heme synthesis and certain cancers)

Abnormal levels of these compounds may reflect mitochondrial stress, cellular damage, inflammation, neoplastic changes, or metabolic disorders. For instance, hematuria (blood in urine) can be confirmed by detecting fluorescence emissions between 450–520 nm, in addition to visible color change.



Hamamatsu's C12666MA mini-spectrometer, fingertip size, ultra-compact spectrometer head integrating MEMS and image sensor technologies

---

## Complementary Sensing-Conductimetry

Spectral methods are enhanced by combining them with conductimetry, which measures the urine's electrical conductivity to assess total ion concentration. This dual approach improves diagnostic accuracy and reduces false negatives, especially when adjusted for sample temperature.

## Advantages of Spectral Methods Over Traditional Techniques

Features	Traditional Methods	Spectral (Optical) Methods
Invasiveness	Non-invasive	Non-invasive
Speed	Minutes to days (e.g., cultures)	Seconds to minutes
Sample preparation	Often required	Minimal or none
Quantitative precision	Limited to lab-based tests	High precision, multi-parametric
Biomarker range	Narrow (specific targets)	Broad (metabolic + structural)
Operator dependency	Moderate to high	Low (automated, software-assisted)
Repeatability	Variable; risk of degradation	High; samples remain intact
Data richness	Low to moderate	High (full spectral fingerprint)



## A New Era in Diagnostics

By integrating visible, infrared, and fluorescence spectrometry with conductimetry into compact analytical platforms, modern urine analyzers are transforming diagnostic medicine. These devices enable clinicians to perform rapid and reliable assessments across a broad spectrum of biomarkers. What was once limited to specialized laboratories can now be conducted in near real-time, offering new possibilities for early disease detection, treatment monitoring, and risk assessment, all without the need for invasive procedures.

As these technologies continue to evolve, the potential for point-of-care and even home-based applications becomes increasingly feasible. Spectral urine analysis is poised to shift diagnostics from the lab bench to the bedside and beyond, making precision healthcare more accessible, efficient, and patient-friendly.

Hamamatsu, with its highly reliable, sensitive, and miniaturized spectral solutions, is playing a pivotal role in this transformation. By enabling more compact and accurate sensing platforms, Hamamatsu is making a significant contribution to the advancement of medical diagnostics, helping to bring innovative, non-invasive tools into everyday clinical practice.

For more information please visit [www.hamamatsu.com](http://www.hamamatsu.com)  
or contact us at [info@hamamatsu.eu](mailto:info@hamamatsu.eu)

---

## References

[1] Senger RS, Kavuru V, Sullivan M, Gouldin A, Lundgren S, Merrifield K, et al. (2019) Spectral characteristics of urine specimens from healthy human volunteers analyzed using Raman chemometric urinalysis (Rametrix). PLoS ONE 14(9): e0222115 [Online]. Available: <https://doi.org/10.1371/journal.pone.0222115>

1. Vigo, F.; Tozzi, A.; Disler, M.; Gisi, A.; Kavadias, V.; Kavadias, T. Vibrational Spectroscopy in Urine Samples as a Medical Tool: Review and Overview on the Current State-of-the-Art. Diagnostics 2023, 13, 27 [Online]. Available: <https://doi.org/10.3390/diagnostics13010027>

2. Li, J.; Wang, C.; Meng, Q.; Hu, Z.; Hu, M.; Zhang, M. MicroRNAs in urine as diagnostic biomarkers for multiple myeloma. Int. J. Lab. Hematol. 2020, 43, 227–234. [CrossRef] [PubMed]

3. Brouwer-Brolsma, E.M.; Brennan, L.; Drevon, C.A.; van Kranen, H.; Manach, C.; Dragsted, L.O.; Roche, H.M.; Andres-Lacueva, C.; Bakker, S.J.L.; Bouwman, J.; et al. Combining traditional dietary assessment methods with novel metabolomics techniques: Present efforts by the Food Biomarker Alliance. Proc. Nutr. Soc. 2017, 76, 619–627. [CrossRef] [PubMed]

4. USense, "Hand-held urine analyzer integrating spectral and electrical measurements," U.S. Patent Application US2024102987A1, Apr. 4, 2024. [Online]. Available: <https://worldwide.espacenet.com/patent/search/family/074141311/publication/US2024102987A1>