

COOLERS FOR PHOTOMULTIPLIER TUBE C10372/C10373 SERIES

The C10372 series and C10373 series are water-cooled thermoelectric coolers designed to reduce thermal electrons emitted from the photocathode of photomultiplier tubes (PMTs) in order to improve signal-to-noise ratio (S/N ratio). The C10372 series further contains an electrostatic and magnetic shield that minimizes the influence of the ambient environment.

The C10372 series are for 28 mm and 51 mm diameter head-on PMTs, while the C10373 series are for MCP-PMTs.

FEATURES

- Thermoelectric cooling using Peltier module
- About -30 °C cooling temperature (with +20 °C cooling water)
- Temperature can be set in the range of -30 °C to 0 °C (dial type)
- Evacuated, double-pane window with heater for frost prevention
- Built-in electrostatic and magnetic shielding (C10372 series)
- Internal protective circuit safeguards Peltier element module in case of low water
- Internal protection circuits prevent output short-circuit, output overvoltage, and excessive temperature rise in power supply
- Stable operation due to a regulated power supply



Left: Power supply Right: Cooled PMT housing



SPECIFICATIONS

COOLED PMT HOUSING

Parameter		Description / Value
Cooling method		Thermoelectric cooling using Peltier module
Heat exchange medium		Water
Amount of cooling water		1 L/min to 3 L/min (water pressure: below 0.3 MPa)
Cooling temperature (with cooling water at +20 °C)		Approx. -30 °C
Temperature controllable range (with cooling water at +20 °C)		-30 °C to 0 °C (continuously adjustable)
Cooling time		Approx. 120 min
Optical window material		Evacuated double-pane synthetic silica window with heater (185 nm to 2200 nm)
Applicable PMTs (sold separately)	C10372 series	28 mm (1-1/8") Dia., 51 mm (2") Dia. head-on
	C10373 series	MCP-PMT (R3809U series, R5916U series) ^(A)
Applicable socket assembly (sold separately)	C10372 series	E2762 series ^(A)
Applicable PMT holder (sold separately)	C10373 series	E3059-500 (R3809U series), E3059-501 (R5916U series)
Operating ambient temperature / Humidity ^(B)		+5 °C to +40 °C / Below 75 %
Storage temperature / Humidity ^{(B)(C)}		-15 °C to +50 °C / Below 80 %
Weight		C10372 5.8 kg / C10373 5.5 kg

NOTE: ^(A)See the section of sold separately ^(B)No condensation ^(C)No cooling water

POWER SUPPLY

Parameter		Description / Value
AC Input voltage	C10372/C10373 ^(D)	AC 100 V (50 Hz / 60 Hz)
	C10372-01/C10373-01 ^(E)	AC 120 V (50 Hz / 60 Hz)
	C10372-02/C10373-02 ^(F)	AC 240 V (50 Hz / 60 Hz)
Power consumption		200 V·A
Output voltage		24 V to 27 V
Output current		4.2 A
Protection circuit		Protective circuits to safeguard Peltier element module in case of low water and to prevent output short-circuit, output overvoltage, and excessive temperature rise in power supply.
Operating ambient temperature / Humidity ^(B)		+5 °C to +40 °C / Below 75 %
Storage temperature / Humidity ^(B)		-15 °C to +50 °C / Below 80 %
Weight		2.1 kg

NOTE: ^(B)No condensation ^(D)For Japan ^(E)For America ^(F)For Europe

COMPONENTS AND ACCESSORIES



- Cooled PMT housing (Including an input window)
- Power supply
- Light shield cap
- Water hose clamps (2 pcs)
- Connection cable (1.5 m)
- AC line cable (2.5 m: -02 type, 2.0 m: other types)
- Socket assembly / PMT holder mounting screw (4 pcs)

* To use these coolers, water hoses with an inner diameter of 15 mm and a water supply line with the matching round faucet are required. Prepare those hoses of the desired length. In addition, prepare a filter for removing impurities such as chlorine ions.

COOLERS FOR PHOTOMULTIPLIER TUBE


C10372/C10373 SERIES

COOLING CHARACTERISTICS




DIMENSIONAL OUTLINE (Unit: mm)

HOUSING




POWER SUPPLY



SOLD SEPARATELY (Unit: mm)


Socket assemblies for C10372 series E2762 series



NOTE: (A)

E2762 series	PMT
E2762-509	R464, R585, R649
E2762-510	R329-02, R331-05, R2257
E2762-511	R316-02, R374, R2228, R5929, R6249

MCP-PMT holder for C10373 series E3059-500/-501



Holder for MCP-PMT	MCP-PMT
E3059-500	R3809U-50 / -52 / -61 / -64
E3059-501	R5916U-50 / -52 / -61

HAMAMATSU PHOTONICS K.K. www.hamamatsu.com

Electron Tube Division

314-5, Shimokanzo, Iwata City, Shizuoka Pref., 438-0193, Japan, Telephone: (81)539-62-5248, Fax: (81)539-62-2205

U.S.A.: HAMAMATSU CORPORATION: 360 Foothill Road, Bridgewater, NJ 08807, U.S.A., Telephone: (1)908-231-0960, Fax: (1)908-231-1218
Germany: HAMAMATSU PHOTONICS DEUTSCHLAND GMBH: Arzbergerstr. 10, 82211 Herrsching am Ammersee, Germany, Telephone: (49)8152-375-0, Fax: (49)8152-265-8 E-mail: info@hamamatsu.de
France: HAMAMATSU PHOTONICS FRANCE S.A.R.L.: 19 Rue du Saule Trapu, Parc du Moulin de Massy, 91882 Massy Cedex, France, Telephone: (33)1 69 53 71 00, Fax: (33)1 69 53 71 10 E-mail: info@hamamatsu.fr
United Kingdom: HAMAMATSU PHOTONICS UK LIMITED: 2 Howard Court, 10 Tewin Road, Welwyn Garden City, Hertfordshire, AL7 1BW, UK, Telephone: (44)1707-294888, Fax: (44)1707-325777 E-mail: info@hamamatsu.co.uk
North Europe: HAMAMATSU PHOTONICS NORDEN AB: Torshamnsgatan 35, 16440 Kista, Sweden, Telephone: (46)8-509-031-00, Fax: (46)8-509-031-01 E-mail: info@hamamatsu.se
Italy: HAMAMATSU PHOTONICS ITALIA S.R.L.: Strada della Moia, 1 int. 6 20044 Arese (Milano), Italy, Telephone: (39)02-93 58 17 33, Fax: (39)02-93 58 17 41 E-mail: info@hamamatsu.it
China: HAMAMATSU PHOTONICS (CHINA) CO., LTD.: 1201, Tower B, Jiaming Center, 27 Dongsanhuang Bellu, Chaoyang District, 100020 Beijing, P.R. China, Telephone: (86)10-6586-6006, Fax: (86)10-6586-2866 E-mail: hpc@hamamatsu.com.cn
Taiwan: HAMAMATSU PHOTONICS TAIWAN CO., LTD.: 13F-1, No.101, Section 2, Gongdao 5th Road, East Dist., Hsinchu City, 300046, Taiwan(R.O.C) Telephone: (886)3-659-0080, Fax: (886)3-659-0081 E-mail: info@hamamatsu.com.tw

TACC1040E04
FEB. 2025 IP