



Flow cells



Flow cells

Hamamatsu provides a wide variety of quartz glass flow cells (cuvettes) designed and fabricated based on our accumulated expertise in glass processing. These flow cells have high dimensional accuracy that ensures the ideal flow rate and flow speed control.

→ P.03



Flow cell assemblies

These flow cell assemblies consist of various component combinations for feeding liquids, and are designed to easily install into analytical instruments such as flow cytometers.

→ P.05




Main application: Flow cytometer

Flow cytometers are a type of cell analyzer that measures the size, number and state of cells. A wide variety of flow cytometers capable of high-precision and high-throughput measurements are available in different sizes and performance specs.

Measurement method

A sample solution containing the target cells is guided to flow through the flow cells while controlling the amount and speed of cell flow so as to align each cell in a single file stream. A laser beam is irradiated onto the flowing cells being aligned, and the resulting generated fluorescence or scattered light is measured.

The flow cells in this way serve as a "pathway" for the target cells and light in the flow cytometer.




High dimensional accuracy ensures precise control of the flow path and optical systems.

Measuring cells with a high degree of accuracy requires precisely controlling flow rate and flow speed of the sample liquid. This requires the flow cells to have high dimensional accuracy since it significantly affects flow control. Hamamatsu flow cells have high dimensional accuracy that ensures precise control of the flow path and also improves flexibility and freedom in designing the optical systems.

Other application

● Blood cell counters



● HPLC



Flow cells

Flow cells J11020 series




Specifications


Parameter	J11020-000-XXX						Unit
Suffix	-23A	-24A	-25A	-33A	-34A	-35A	
Flow channel size	0.150 × 0.150	0.200 × 0.200	0.250 × 0.250	0.150 × 0.150	0.200 × 0.200	0.250 × 0.250	mm
Optical effective areas	5.0 × 2.0			Φ2.0 / 2.0 × 0.6			mm
Material	Synthetic silica						—
Operating ambient temperature	+15 to +35						°C
Operating ambient relative humidity	<85						%RH
Storage temperature	0 to +45						°C
Weight	0.87			0.57			g

Dimensional outline examples (Unit: mm)

- Standard type (J11020-000-25A)




- Thin cover glass type (J11020-000-35A)




* The external dimensions of the J11020-000-23A/-24A/-33A/-34A are the same as the -25A and -35A except for the channel size. See the specification table to find their channel size.

Customization examples


We welcome requests for different shapes and custom options.
Customization tasks not listed on this page will also be available, so please contact us for details.




Shape custom: with step




Shape custom: with lens




Shape custom: thin cover glass




Shape custom: Large channel




Other custom: Marking




Other custom: With grinding (Pattern1)



Other custom: With grinding (Pattern2)




Material custom: Synthetic black quartz




Other custom: AR coating

● Transmittance of synthetic black quartz



* The values on the above chart are sample values only
* Cannot be used for the UV spectral range.


● Optical thin film coating



Ex.: AR (Anti-Reflection) coating: reflectance characteristics
* The values on the above chart are sample values only

Flow cell assemblies

Flow cell assemblies J12800 series




Specifications

Parameter		J12800-000-XXX						Unit	
Suffix		-A3A	-A4A	-A5A	-B3A	-B4A	-B5A		
Flow cell	Flow channel size	0.15 × 0.15	0.20 × 0.20	0.25 × 0.25	0.15 × 0.15	0.20 × 0.20	0.25 × 0.25	mm	
	Optical effective areas	5.0 × 2.0			Φ2.0			mm	
	Material	Synthetic silica						—	
	Pressure resistance	Max.	450						kPa
Assembly	Recommended input sheath flow rate	7.0	13.5	19.0	7.0	13.5	19.0	ml/min	
	Input sheath flow rate	Max.	9.5	15	24	9.5	15	24	ml/min
	Operating sheath pressure		10 to 180	10 to 120	10 to 100	10 to 180	10 to 120	10 to 100	kPa
	Pressure resistance	Max.	450						kPa
	Connector ①	Sample inlet port	IDEX P-646						—
		Sheath inlet port	IDEX P-647						—
		Waste port	IDEX P-646						—
		Reverse flow port	IDEX P-646						—
	Sample needle inside diameter		0.15						mm
	Operating ambient temperature		+15 to +35						°C
Operating ambient relative humidity ②		<85						%RH	
Storage temperature		0 to +45						°C	
Weight		60.6			60.3			g	

① See P.6 for applicable tubes. ② No condensation.


Dimensional outline (Unit: mm)

• Standard type (J12800-000-A5A)



* The external dimensions of the J12800-000-A3A/-A4A are the same as the -A5A except for the channel size. See the specification table to find their channel size.

● Thin cover glass type (J12800-000-B5A)



* The external dimensions of the J12800-000-B3A/-B4A are the same as the -B5A except for the channel size. See the specification table to find their channel size.

● Tube size and fitting screw list

Item	Tube size	Fitting screw
Sheath inlet	3.20 mm (1/8")	1/4-28UNF
Bubble release	1.58 mm (1/16")	1/4-28UNF
Waste	1.58 mm (1/16")	1/4-28UNF
Sample inlet	1.58 mm (1/16")	1/4-28UNF

Q&A

Flow cells questions

Q.01 Is it possible to change the flow cell specs such as flow channel size and wall thickness?

A.01 Yes, it is possible to customize the specs. Please contact us to discuss your specifications, applications, and requirements.

Q.02 I need a 3D model for optical design. Can you provide it?

A.02 Yes, we can provide it.

Q.03 Can you give us information on the basic optical properties of the flow cell glass such as transmittance?

A.03 Yes, that information is available.

Q.04 I'd like to consult with you about a customization need. Do you require any information for the consultation?

A.04 It would make things easier to discuss if you could provide us with a simple drawing or figure for the customization.

Q.05 Are flow cells made from material other than quartz glass available?

A.05 We are sorry but these are not available in other than quartz glass.

Flow cell assemblies questions

Q.01 Can you fabricate a flow cell assembly using a flow cell with an AR lens or by processing black quartz?

A.01 It may be possible in some cases. Please contact us for details.

Q.02 Is it possible to drill mounting holes in the plastic or metal parts of a flow cell assembly?

A.02 Yes, it is possible. Please consult with us about this.

Q.03 Can I replace or remove an already assembled flow cell?

A.03 No, you can't because the flow cell is assembled by using adhesives.

Subject to local technical requirements and regulations, availability of products included in this promotional material may vary. Please consult with our sales office.
Information furnished by HAMAMATSU is believed to be reliable. However, no responsibility is assumed for possible inaccuracies or omissions. Specifications are subject to change without notice. No patent rights are granted to any of the circuits described herein. ©2023 Hamamatsu Photonics K.K.

HAMAMATSU PHOTONICS K.K. www.hamamatsu.com

Electron Tube Division

314-5, Shimokanzo, Iwata City, Shizuoka Pref., 438-0193, Japan, Telephone: (81)539/62-5248, Fax: (81)539/62-2205

U.S.A.: HAMAMATSU CORPORATION: 360 Foothill Road, Bridgewater, NJ 08807, U.S.A., Telephone: (1)908-231-0960, Fax: (1)908-231-1218 E-mail: usa@hamamatsu.com

Germany: HAMAMATSU PHOTONICS DEUTSCHLAND GMBH.: Arzbergerstr. 10, 82211 Hersching am Ammersee, Germany, Telephone: (49)8152-375-0, Fax: (49)8152-265-8 E-mail: info@hamamatsu.de

France: HAMAMATSU PHOTONICS FRANCE S.A.R.L.: 19 Rue du Saule Trapu, Parc du Moulin de Massy, 91882 Massy Cedex, France, Telephone: (33)1 69 53 71 00, Fax: (33)1 69 53 71 10 E-mail: infos@hamamatsu.fr

United Kingdom: HAMAMATSU PHOTONICS UK LIMITED: 2 Howard Court, 10 Tewin Road, Welwyn Garden City, Hertfordshire, AL7 1BW, UK, Telephone: (44)1707-294888, Fax: (44)1707-325777 E-mail: info@hamamatsu.co.uk

North Europe: HAMAMATSU PHOTONICS NORDEN AB: Torshamnsgatan 35, 16440 Kista, Sweden, Telephone: (46)8-509-031-00, Fax: (46)8-509-031-01 E-mail: info@hamamatsu.se

Italy: HAMAMATSU PHOTONICS ITALIA S.R.L.: Strada della Moia, 1 int. 6 20044 Arese (Milano), Italy, Telephone: (39)02-93 58 17 33, Fax: (39)02-93 58 17 41 E-mail: info@hamamatsu.it

China: HAMAMATSU PHOTONICS (CHINA) CO., LTD.: 1201, Tower B, Jiaxing Center, 27 Dongsanhuan Beilu, Chaoyang District, 100020 Beijing, P.R. China, Telephone: (86)10-6586-6006, Fax: (86)10-6586-2866 E-mail: hpc@hamamatsu.com.cn

Taiwan: HAMAMATSU PHOTONICS TAIWAN CO., LTD.: 8F-3, No.158, Section 2, Gongdao 5th Road, East District, Hsinchu, 300, Taiwan R.O.C. Telephone: (886)3-659-0080, Fax: (886)3-659-0081 E-mail: info@hamamatsu.com.tw

TPMZ1011E05
MAY 2023