

FEATURES

- High quantum efficiency: GaAsP / GaAs photocathode
- Built-in thermoelectric cooler
- Low noise
- For photon counting



▲H16721 + A7423 (Heatsink with Fan)

SPECIFICATIONS

(at +25 °C)

Parameter		H16721-40	H16721-50	Unit
Input voltage		+4.5 to +5.5		V
Max. input voltage for main unit		+6		V
Max. input current for main unit		50		mA
Max. input voltage for thermoelectric cooler		2.6		V
Max. input current for thermoelectric cooler		2.2		A
Effective area		φ5		mm
Photocathode material		GaAsP	GaAs	—
Spectral response		300 to 740	380 to 900	nm
Peak count sensitivity wavelength		610	730	nm
Count sensitivity	400 nm	5.3×10^5	—	s ⁻¹ ·pW ⁻¹
	500 nm	9.0×10^5	2.6×10^5	
	600 nm	10.3×10^5	4.5×10^5	
	700 nm	8.5×10^5	5.1×10^5	
	800 nm	—	4.8×10^5	
Count linearity ^①		1.5×10^6		s ⁻¹
Dark count ^{②③}	Typ.	100	125	s ⁻¹
	Max.	300	375	
Pulse-pair resolution		70		ns
Output pulse width		30		ns
Output pulse height ^④	Min.	3.0		V
	Typ.	3.6		
Recommended load resistance		50		Ω
Signal output logic		Positive logic		—
Operating ambient temperature ^⑤		+5 to +35		°C
Storage temperature ^⑤		-20 to +50		°C
Weight		Approx. 395		g

①Random pulse, at 10 % count loss

②PMT setting temperature 0 °C, used with C8137 and A7423

③After 30 min storage in darkness

④With input voltage +5 V, Load resistance 50 Ω and coaxial cable RG-174/U (450 mm)

⑤No condensation

COOLING SPECIFICATIONS

Parameter	H16721-40 / H16721-50	Unit
Cooling method	Thermoelectric cooling	—
Max. cooling temperature (ΔT) ^⑥	35	°C
Time to reach Max. cooling temperature ^⑥	Approx. 5	min

⑥Input current to thermoelectric cooler = 2 A

PHOTON COUNTING HEAD H16721 SERIES

Figure 1: Count sensitivity




Figure 2: Dark count




Figure 3: Block diagram






Figure 4: Dimensional outline (Unit: mm)



OPTIONS FOR H16721 SERIES



● Heatsink with fan A7423

The temperature of the H16721 outer case rises due to the thermoelectric cooler housed in the case. The A7423 heatsink efficiently radiates away this heat to prevent a temperature rise in the H16721. The A7423 can be easily installed onto the H16721 with four M3 screws. Apply a coat of heat conductive grease onto the joint surface shared by the H16721 and A7423.

Parameter	Value	Unit	
Input voltage	12	V	
Input current	During lock	140	mA
	During operation	90	mA
Operating voltage	10.2 to 13.8	V	
Weight	128	g	

● Signal cable E1168-05

This signal cable comes attached to a BNC connector for easily connecting the H16721 to external equipment.

● Optical fiber adapter (FC type) A7412

The A7412 is an FC type optical fiber connector that attaches to the light input window of the H16721. The A7412 can easily be secured in place with four M2 screws.

● C-mount adapter A7413

The A7413 mount adapter is used when a C-mount lens protruding 4 mm or more from the flange-back must be installed onto the H16721.

● Power supply unit with temperature control C8137


The C8137 is a power supply unit with a temperature control function. Just connecting to an AC source of 100 V to 240 V generates the output voltages for the thermoelectric cooler and the A7423 fan, needed for operating the H16721. The photomultiplier tube temperature can be maintained to 0 °C by monitoring the thermistor and regulating the output current for the thermoelectric cooler.

Parameter	Value	Unit
Max. cooling temperature (ΔT)	35	°C
Setting cooling temperature (preset at factory)	0	°C
AC input voltage	100 to 240	V
Input voltage frequency	50 / 60	Hz
Power consumption	30	V·A
Main circuit output voltage	+5	V
Max. current for thermoelectric cooler	2.2	A
Output voltage for fan	12	V
Weight	1.1	kg


PHOTON COUNTING HEAD H16721 SERIES

OPTIONS (Unit: mm)

●Heatsink with fan A7423




TOP VIEW



SIDE VIEW

TACCA0188ED

●Optical fiber adapter (FC type) A7412




FRONT VIEW

SIDE VIEW


TACCA0190EB

●C-mount adapter A7413



TACCA0191EA




●Power supply unit with temperature control C8137



FRONT VIEW


SIDE VIEW

REAR VIEW



TACCA0237EB

●Signal cable E1168-05



TACCA0148EA

HAMAMATSU PHOTONICS K.K. www.hamamatsu.com

Electron Tube Division

314-5, Shimokanzo, Iwata City, Shizuoka Pref., 438-0193, Japan, Telephone: (81)539/62-5248, Fax: (81)539/62-2205

U.S.A.: HAMAMATSU CORPORATION: 360 Foothill Road, Bridgewater, NJ 08807, U.S.A., Telephone: (1)908-231-0960, Fax: (1)908-231-1218 E-mail: usa@hamamatsu.com

Germany: HAMAMATSU PHOTONICS DEUTSCHLAND GMBH.: Arzbergerstr. 10, 82211 Herrsching am Ammersee, Germany, Telephone: (49)8152-375-0, Fax: (49)8152-265-8 E-mail: info@hamamatsu.de

France: HAMAMATSU PHOTONICS FRANCE S.A.R.L.: 19 Rue du Saule Trapu, Parc du Moulin de Massy, 91882 Massy Cedex, France, Telephone: (33)1 69 53 71 10 E-mail: info@hamamatsu.fr

United Kingdom: HAMAMATSU PHOTONICS UK LIMITED: 2 Howard Court, 10 Tewin Road, Welwyn Garden City, Hertfordshire, AL7 1BW, UK, Telephone: (44)1707-294888, Fax: (44)1707-325777 E-mail: info@hamamatsu.co.uk

North Europe: HAMAMATSU PHOTONICS NORDEN AB: Torshamnsgatan 35 16440 Kista, Sweden, Telephone: (46)8-509 031 00, Fax: (46)8-509 031 01 E-mail: info@hamamatsu.se

Italy: HAMAMATSU PHOTONICS ITALIA S.R.L.: Strada della Moia, 1 int. 6, 20044 Arese (Milano), Italy, Telephone: (39)02-93 58 17 33, Fax: (39)02-93 58 17 41 E-mail: info@hamamatsu.it

China: HAMAMATSU PHOTONICS (CHINA) CO., LTD.: 1201 Tower B, Jiaming Center, 27 Dongsanhuan Bellu, Chaoyang District, 100020 Beijing, P.R. China, Telephone: (86)10-6586-6006, Fax: (86)10-6586-2866 E-mail: hpc@hamamatsu.com.cn

Taiwan: HAMAMATSU PHOTONICS TAIWAN CO., LTD.: 8F-3, No.158, Section 2, Gongdao 5th Road, East District, Hsinchu, 300, Taiwan R.O.C. Telephone: (886)3-659-0080, Fax: (886)3-659-0081 E-mail: info@hamamatsu.com.tw

TPMO1121E01
OCT. 2022 IP