

FEATURES

- Wide spectral response
 R3788 185 nm to 750 nm
 R4332 160 nm to 750 nm
- High cathode sensitivity
 Luminous 120 μ A/lm Typ.
 Radiant at 420 nm 90 mA/W Typ.
 Quantum efficiency at 210 nm 40 % Typ. (R4332)
- High anode sensitivity (at 1000 V)
 Luminous 1200 A/lm Typ.
 Radiant at 420 nm 9.0×10^5 A/W Typ.



▲Left: R3788, Right: R4332

APPLICATIONS


- Fluorescence spectrophotometers
- Emission spectrophotometers
- Atomic absorption spectrophotometers

SPECIFICATIONS

GENERAL

Parameter	Description / Value	Unit
Spectral response	R3788	185 to 750 nm
	R4332	160 to 750 nm
Wavelength of maximum response	420	nm
Photocathode	Material	Bialkali
	Minimum effective area	8 × 24
Window material	R3788	UV glass
	R4332	Silica glass
Dynode	Secondary emitting surface	Bialkali
	Structure	Circular-cage
	Number of stages	9
Direct Interelectrode capacitances	Anode to last dynode	4 pF
	Anode to all other electrodes	6 pF
Base	11-pin base JEDEC No. B11-88	—
Weight	Approx. 45	g
Operating ambient temperature	-30 to +50	°C
Storage temperature	-30 to +50	°C
Suitable socket	E678-11A (Sold separately)	—
Suitable socket assembly	E717-63 (Sold separately)	—
	E717-74 (Sold separately)	—

Figure 1: Typical spectral response



PHOTOMULTIPLIER TUBES R3788, R4332

MAXIMUM RATINGS (Absolute maximum values)

Parameter		Value	Unit
Supply voltage	Between anode and cathode	1250	V
	Between anode and last dynode	250	V
Average anode current ^A		0.1	mA

CHARACTERISTICS (at 25 °C)

Parameter		R3788			R4332			Unit	
		Min.	Typ.	Max.	Min.	Typ.	Max.		
Cathode sensitivity	Quantum efficiency (at peak wavelength)	—	30 (at 250 nm)	—	—	40 (at 210 nm)	—	%	
	Luminous ^B	100	120	—	100	120	—	μA/lm	
	Radiant	at 194 nm	—	31	—	—	60	—	mA/W
		at 210 nm	—	50	—	—	68	—	mA/W
		at 420 nm	—	90	—	—	90	—	mA/W
	Red / White ratio ^C	—	0.01	—	—	0.01	—	—	
Blue sensitivity index ^D	—	10	—	—	10	—	—		
Anode sensitivity	Luminous ^E	500	1200	—	500	1200	—	A/lm	
	Radiant	at 194 nm	—	3.1 × 10 ⁵	—	—	6.0 × 10 ⁵	—	A/W
		at 210 nm	—	5.0 × 10 ⁵	—	—	6.8 × 10 ⁵	—	A/W
		at 420 nm	—	9.0 × 10 ⁵	—	—	9.0 × 10 ⁵	—	A/W
Gain ^E		—	1.0 × 10 ⁷	—	—	1.0 × 10 ⁷	—	—	
Anode dark current ^F (After 30 min storage in darkness)		—	5	50	—	5	50	nA	
ENI (Equivalent Noise Input) ^G		—	1.4 × 10 ⁻¹⁶	—	—	1.4 × 10 ⁻¹⁶	—	W	
Time response ^E	Anode pulse rise time ^H	—	2.2	—	—	2.2	—	ns	
	Electron transit time ^I	—	22	—	—	22	—	ns	
	Transit time spread (T.T.S.) ^J	—	1.2	—	—	1.2	—	ns	
Anode current stability ^K	Light hysteresis	—	0.1	—	—	0.1	—	%	
	Voltage hysteresis	—	1.0	—	—	1.0	—	%	

NOTES

- A: Averaged over any interval of 30 s maximum.
 B: The light source is a tungsten filament lamp operated at a distribution temperature of 2856 K. Supply voltage is 100 V between the cathode and all other electrodes connected together as anode.
 C: Red / White ratio is the quotient of the cathode current measured using a red filter interposed between the light source and the tube by the cathode current measured with the filter removed under the same conditions as Note B.
 D: The value is cathode output current when a blue filter is interposed between the light source and the tube under the same condition as Note B.
 E: Measured with the same light source as Note B and with the voltage distribution ratio shown in Table 1 below.

Table 1: Voltage distribution ratio

Electrodes	K	Dy1	Dy2	Dy3	Dy4	Dy5	Dy6	Dy7	Dy8	Dy9	P
Distribution ratio	1	1	1	1	1	1	1	1	1	1	1

Supply voltage: 1000 V, K: Cathode, Dy: Dynode, P: Anode


- F: Measured with the same supply voltage and voltage distribution ratio as Note E after removal of light.
 G: ENI is an indication of the photon-limited signal-to-noise ratio. It refers to the amount of light in watts to produce a signal-to-noise ratio of unity in the output of a photomultiplier tube.

$$ENI = \frac{\sqrt{2q \cdot I_{db} \cdot G \cdot \Delta f}}{S}$$

- where q = Electronic charge (1.60 × 10⁻¹⁹ coulomb).
 I_{db} = Anode dark current(after 30 min storage) in amperes.
 G = Gain.
 Δf = Bandwidth of the system in hertz. 1 hertz is used.
 S = Anode radiant sensitivity in A/W at the wavelength of peak response.

- H: The rise time is the time for the output pulse to rise from 10 % to 90 % of the peak amplitude when the entire photocathode is illuminated by a delta function light pulse.
 I: The electron transit time is the interval between the arrival of delta function light pulse at the entrance window of the tube and the time when the anode output reaches the peak amplitude. In measurement, the whole photocathode is illuminated.
 J: Also called transit time jitter. This is the fluctuation in electron transit time between individual pulses in the signal photoelectron mode, and may be defined as the FWHM of the frequency distribution of electron transit times.
 K: Hysteresis is temporary instability in anode current after light and voltage are applied.

$$\text{Hysteresis} = \frac{I_{\max.} - I_{\min.}}{I_i} \times 100 (\%)$$



TPMSB0002EA

(1) Light hysteresis

The tube is operated at 750 V with an anode current of 1 μA for 5 minutes. The light is then removed from the tube for a minute. The tube is then re-illuminated by the previous light level for a minute to measure the variation.

(2) Voltage hysteresis

The tube is operated at 300 V with an anode current of 0.1 μA for 5 minutes. The light is then removed from the tube and the supply voltage is quickly increased to 800 V. After a minute, the supply voltage is then reduced to the previous value and the tube is re-illuminated for a minute to measure the variation.

Figure 2: Anode luminous sensitivity and gain characteristics




Figure 3: Typical time response




Figure 4: Typical ENI with wavelength





Figure 5: Typical EADCI (Equivalent Anode Dark Current Input) vs. supply voltage



PHOTOMULTIPLIER TUBES R3788, R4332

Figure 6: Dimensional outline and basing diagram (Unit: mm)




Figure 7: Socket (Unit: mm) Sold separately






Figure 8: D type socket assembly (Unit: mm) Sold separately

E717-63



E717-74



* Hamamatsu also provides C13890 series compact high voltage power supplies and C6270 series DP type socket assemblies which incorporate a DC to DC converter type high voltage power supply.

Warning—Personal Safety Hazards

Electrical Shock—Operating voltages applied to this device present a shock hazard.

HAMAMATSU PHOTONICS K.K. www.hamamatsu.com

Electron Tube Division

314-5, Shimokanzo, Iwata City, Shizuoka Pref., 438-0193, Japan, Telephone: (81)539/62-5248, Fax: (81)539/62-2205

U.S.A.: Hamamatsu Corporation, 360 Foothill Road, Bridgewater, NJ 08807, U.S.A., Telephone: (1)908-231-0960, Fax: (1)908-231-1218 E-mail: usa@hamamatsu.com

Germany: Hamamatsu Photonics Deutschland GmbH, Arzbergerstr. 10, D-82211 Herrsching am Ammersee, Germany, Telephone: (49)8152-375-0, Fax: (49)8152-265-8 E-mail: info@hamamatsu.de

France: Hamamatsu Photonics France S.A.R.L., 19, Rue du Saule Trapu, Parc du Moulin de Massy, 91882 Massy Cedex, France, Telephone: (33)1 69 53 71 00, Fax: (33)1 69 53 71 10 E-mail: infos@hamamatsu.fr

United Kingdom: Hamamatsu Photonics UK Limited, 2 Howard Court, 10 Tewin Road, Welwyn Garden City, Hertfordshire AL7 1BW, UK, Telephone: (44)1707-294888, Fax: (44)1707-325777 E-mail: info@hamamatsu.co.uk

North Europe: Hamamatsu Photonics Norden AB, Torshamnsgatan 35 16440 Kista, Sweden, Telephone: (46)8-509 031 00, Fax: (46)8-509 031 01 E-mail: info@hamamatsu.se

Italy: Hamamatsu Photonics Italia S.r.l., Strada della Moia, 1 int. 6, 20044 Arese (Milano), Italy, Telephone: (39)02-93 58 17 33, Fax: (39)02-93 58 17 41 E-mail: info@hamamatsu.it

China: Hamamatsu Photonics (China) Co., Ltd., 1201 Tower B, Jiaming Center, 27 Dongshanhuan Bellu, Chaoyang District, 100020 Beijing, P.R. China, Telephone: (86)10-6586-6006, Fax: (86)10-6586-2866 E-mail: hpc@hamamatsu.com.cn

Taiwan: Hamamatsu Photonics Taiwan Co., Ltd., 8F-3, No.158, Section 2, Gongdao 5th Road, East District, Hsinchu, 300, Taiwan R.O.C. Telephone: (886)3-659-0080, Fax: (886)3-659-0081 E-mail: info@hamamatsu.com.tw

TPMS1021E03
JUN. 2021 IP