

For vacuum ultraviolet light detection
Cs-Te (R8486), Cs-I (R8487) photocathode, MgF₂ window,
28 mm (1-1/8 inch) diameter, 9-stage, Side-on type

FEATURES

- Sensitivity in the vacuum ultraviolet region
 - R8486.....115 to 320 nm
 - R8487.....115 to 195 nm
- High quantum efficiency (at 121.6 nm)
 - R8486.....22.5 % (Typ.)
 - R8487.....26.0 % (Typ.)
- High anode sensitivity
 - R8486 (at 254 nm)..... 5.2×10^5 A/W (Typ.)
 - R8487 (at 121.6 nm)..... 1.0×10^5 A/W (Typ.)



APPLICATIONS


- Emission spectroscopy, etc.

SPECIFICATIONS

GENERAL

Parameter		R8486	R8487	Unit
Spectral response		115 to 320	115 to 195	nm
Wavelength of maximum response		230	130	nm
Photocathode material		Cs-Te	Cs-I	—
Window material		MgF ₂		—
Minimum effective area		8 × 12		mm
Dynode	Structure	Circular-cage		—
	Number of stage	9		—
	Material	Sb-Cs		—
Direct interelectrode capacitances	Anode to dynode No.9	Approx. 4		pF
	Anode to all other electrodes	Approx. 6		pF
Base		11-pin base JEDEC No. B11-88		—
Weight		48.5		g
Suitable socket		E678-11A (sold separately)		—
Operating ambient temperature		-30 to +50		°C
Storage temperature		-30 to +50		°C

Figure 1: Typical spectral response



PHOTOMULTIPLIER TUBES R8486, R8487

MAXIMUM RATINGS (Absolute maximum values)

Parameter	Rating	Unit
Supply voltage		
Between anode and cathode	1250	V
Between anode and last dynode	250	V
Average anode current ^(A)	0.1	mA

CHARACTERISTICS (at 25 °C)

Parameter	R8486	R8487	Unit
Cathode sensitivity			
Quantum efficiency at 121.6 nm	22.5	26.0	%
at 254 nm	25.0	—	%
Anode sensitivity ^(B)			
Radiant at 121.6 nm	—	1.0 × 10 ⁵	A/W
at 254 nm	5.2 × 10 ⁵	—	A/W
Gain	1.0 × 10 ⁷	3.9 × 10 ⁶	—
Anode dark current (After 30 minute storage in darkness) ^(B)	1.0	0.1	nA
ENI (Equivalent Noise Input) ^(C) at 121.6 nm	—	1.12 × 10 ⁻¹⁶	W
at 254 nm	1.09 × 10 ⁻¹⁶	—	W
Time response			
Anode pulse rise time ^(D)	2.2	2.2	ns
Electron transit time ^(E)	22	22	ns
Transit time spread ^(F)	1.2	1.2	ns

NOTES

- (A): Averaged over any interval of 30 seconds maximum.
 (B): Measured with the voltage distribution ratio shown in Table 1 below.

Table 1: Voltage distribution ratio

Electrode	K	Dy1	Dy2	Dy3	Dy4	Dy5	Dy6	Dy7	Dy8	Dy9	P
Distribution ratio	1	1	1	1	1	1	1	1	1	1	1

Supply voltage=1000 V

K: Cathode Dy: Dynode P: Anode

- (C): ENI is an indication of the photon-limited signal-to-noise ratio. It refers to the amount of light in watts to produce a signal-to-noise ratio of unity in the output of a photomultiplier tube.

$$ENI = \frac{\sqrt{2q \cdot I_{db} \cdot g \cdot \Delta f}}{S}$$

- where q = Electronic charge (1.60 × 10⁻¹⁹ coulomb).
 I_{db} = Anode dark current in amperes after 30 minutes storage in darkness.
 g = Gain.
 Δf = Bandwidth of the system in hertz.
 S = Anode radiant sensitivity in amperes per watt at the wavelength of interest.

- (D): The rise time is the time for the output pulse to rise from 10 % to 90 % of the peak amplitude when the entire photocathode is illuminated by a delta function light pulse.
 (E): The electron transit time is the interval between the arrival of delta function light pulse at the entrance window of the tube and the time when the anode output reaches the peak amplitude. In measurement, the whole photocathode is illuminated.
 (F): Also called transit time jitter. This is the fluctuation in electron transit time between individual pulses in the signal photoelectron mode, and may be defined as the FWHM of the frequency distribution of electron transit times.

Figure 2: Typical gain and anode radiant sensitivity




Figure 3: Typical time response




Figure 4: Dimensional outline and basing diagram (Unit: mm)





Figure 5: Socket E678-11A (Sold separately) (Unit: mm)



NOTE: There is a 2 mm diameter hole to exhaust inner air on the plastic base.

PHOTOMULTIPLIER TUBES R8486, R8487

Warning—Personal safety hazards

Electrical shock—Operating voltages applies to this device present a shock hazard.

HAMAMATSU PHOTONICS K.K. www.hamamatsu.com

Electron Tube Division
314-5, Shimokanzo, Iwata City, Shizuoka Pref., 438-0193, Japan, Telephone: (81)539/62-5248, Fax: (81)539/62-2205
U.S.A.: HAMAMATSU CORPORATION: 360 Foothill Road, Bridgewater, NJ 08807, U.S.A., Telephone: (1)908-231-0960, Fax: (1)908-231-1218
Germany: HAMAMATSU PHOTONICS DEUTSCHLAND GMBH.: Arzbergerstr. 10, 82211 Herrsching am Ammersee, Germany, Telephone: (49)8152-375-0, Fax: (49)8152-265-8 E-mail: info@hamamatsu.de
France: HAMAMATSU PHOTONICS FRANCE S.A.R.L.: 19 Rue du Saule Trapu, Parc du Moulin de Massy, 91882 Massy Cedex, France, Telephone: (33)1 69 53 71 00, Fax: (33)1 69 53 71 10 E-mail: infos@hamamatsu.fr
United Kingdom: HAMAMATSU PHOTONICS UK LIMITED: 2 Howard Court, 10 Tewin Road, Welwyn Garden City, Hertfordshire, AL7 1BW, UK, Telephone: (44)1707-294888, Fax: (44)1707-325777 E-mail: info@hamamatsu.co.uk
North Europe: HAMAMATSU PHOTONICS NORDEN AB: Torshamnsgatan 35, 16440 Kista, Sweden, Telephone: (46)8-509-031-00, Fax: (46)8-509-031-01 E-mail: info@hamamatsu.se
Italy: HAMAMATSU PHOTONICS ITALIA S.R.L.: Strada della Moia, 1 int. 6 20044 Arese (Milano), Italy, Telephone: (39)02-93 58 17 33, Fax: (39)02-93 58 17 41 E-mail: info@hamamatsu.it
China: HAMAMATSU PHOTONICS (CHINA) CO., LTD.: 1201, Tower B, Jiaming Center, 27 Dongsanhuan Beilu, Chaoyang District, 100020 Beijing, P.R. China, Telephone: (86)10-6586-6006, Fax: (86)10-6586-2866 E-mail: hpc@hamamatsu.com.cn
Taiwan: HAMAMATSU PHOTONICS TAIWAN CO., LTD.: 13F-1, No.101, Section 2, Gongdao 5th Road, East Dist., Hsinchu City, 300046, Taiwan(R.O.C) Telephone: (886)3-659-0080, Fax: (886)3-659-0081 E-mail: info@hamamatsu.com.tw

TPMS1070E05
SEP. 2024 IP