

Lateral flow readers THETA[®]



High sensitivity and high reproducibility with ease of use

Hamamatsu Photonics supplies lateral flow readers for R&D work and quality control of immunochromatographic reagents as well as for OEM applications. Lateral flow readers make rapid quantitative measurements of color and fluorescence intensity of immunochromatographic reagents using red/blue-based color particles or fluorescent particles as labels. Our lateral flow readers utilize the latest optical measurement technologies where we never stop exploring and pushing their limits. Try out our lateral flow readers for yourself and see the performance we can offer as a leading manufacturer of optical sensors!

Reagent development and reagent quality control



We offer powerful tools for development and quality control of immunochromatographic reagents. Our product lineup includes a variety of models to ideally match different detection targets and sensitivities, so you can select the optimal model for the target reagent you want to develop and use.

→ P.3 ~

Supports development of lateral flow readers



We supply a 'core engine' for lateral flow readers as a single separate unit utilizable for developing lateral flow readers that deliver high-sensitivity line detection performance.

→ P.5

OEM supply of lateral flow readers



We design and develop lateral flow readers optimized for your specific reagents. Based on the C16723, we will work to mass produce high-quality custom readers and also design specialized devices including interfaces.



→ P.7

Features

Based on optical measurement technologies we have carefully built up over many years, we are continually improving our lateral flow readers to enhance their core performance such as sensitivity and reproducibility. We also weave a variety of underlying technologies into that work to develop and manufacture sophisticated devices having even more highly optimized designs.

● High sensitivity


The C16723 series lateral flow readers detect test lines that are virtually impossible for human eyes to identify. The C16723-60 in particular is designed to support the europium chemical element, providing higher sensitivity than other C16723 series models.



* Measurement sample: standard sample

● High reproducibility

The C16723 series suppresses measurement variations to lower than 3%. This allows accurate evaluation of lot-to-lot differences in reagent performance and measurement variations over time.



CV: Coefficient of variation indicating analysis accuracy (ratio of the standard deviation to the mean during repeated measurements)

● Reagent compatibility

We have vast experience in detecting a wide range of labeled samples using both color-developing and fluorescent reagents. The C16723 series utilizes a flexible universal holder that accommodates various types of reagent trays.



Performance comparison of measurement methods


Item	Visual inspection	Ordinary reader	lateral flow reader
Sensitivity	Low Judgment results differ depending on the person.	Medium Provides higher sensitivity than visual inspections	High Can detect extremely faint color and fluorescence
Reproducibility	No No reproducibility, since judgment cannot be quantified.	Low Judgment results may vary depending on the reader's accuracy.	High High reproducibility with CV (coefficient of variation) kept within 3%
Reagent compatibility	Limited Only reagents that can be checked visually.	Limited Measurable reagents are limited.	Wide compatibility Provides a lineup of readers compatible with different reagents and various shapes.

Functions

The software for the lateral flow reader series has two measurement modes that are selectable to match your task. Waveforms of measurement results are automatically acquired and their peaks then analyzed and displayed based on those waveforms.

* A PC is required to run the dedicated software. Customers must prepare a PC as it is not supplied with our lateral flow readers.


● Example of automatically acquired waveforms



Lot measurement mode

This mode continuously measures multiple reagents of the same lot (same type). Measurement results are saved in CSV format and can be viewed in their order of measurement

● Continuous measurement example



● Acquired analysis data example (CSV)

Sample No.	Peak-①		Peak-②	
	PeakX	mABS	PeakX	mABS
1	292	15.4	699	199.5
2	298	15.1	706	199.0
3	297	15.2	705	198.9
4	296	15.1	704	199.1
5	300	15.2	709	199.0

Time-course measurement mode

This mode measures the reaction process of one reagent and follows its changes over time. This is effective in checking the reaction speed of each reagent and in analyzing time-course changes in color development and fluorescence intensities.

● Time course measurement example



● Acquired analysis data example (time-course changes in peak absorbance)



Advantages over previous units

- Rapid measurement time: about one-half that of the previous units
- Supports a broader range of samples sizes
- Pass-fail check of color samples
- Display screen (for easy verification of measurement count, etc.)
- Improved versatility and operability via upgrades of measurement software

C16723 series




The C16723 series is a family of lateral flow readers optimized for reagent development and quality control. The lineup includes 5 models: C16723-10 and -11 that utilize the absorption method; and C16723-50, -51 and -60 that utilize the fluorescence method. All models have excellent sensitivity, reproducibility and reliability.

Item	Description / value				
Suffix	-10	-11	-50	-51	-60
Target label	Red-based color lines	Blue-based color lines	Europium	FITC	Europium
Measurement target ①	Color labels		Fluorescent labels		
Measurement method	Absorption method		Fluorescence method		High-sensitivity Fluorescence method
Light source	Green LED	Red LED	UV LED	Blue LED	UV LED

①Two colors cannot be measured simultaneously.


Waveform analysis

Absorption profile waveform image (C16723-10/-11)




Color development intensity is automatically calculated by dedicated software.

Fluorescence profile waveform image (C16723-50/-51/-60)



Fluorescence intensity is automatically calculated by dedicated software.

Dimensional outlines (unit: mm)



Specifications

Main unit

Item	C16723-10	C16723-11	C16723-50	C16723-51	C16723-60	Unit
Target label	Red-based color lines	Blue-based color lines	Europium	FITC	Europium	-
Features	Absorption method		Fluorescence method		High-sensitivity Fluorescence method	-
Power requirements	Max.	5				V
Interface	USB2.0					-
Light source	Green LED	Red LED	UV LED	Blue LED	UV LED	-
Detector	Silicon photodiode					-
Dimensions (W × H × D) ^①	100 × 100 × 225					mm

① Not including protruding parts.

Software (for PC operation)

Item	Absorption method type	Fluorescence method type	Unit
Supported OS	Windows® 10 (32 bit / 64 bit), Windows® 11		-
Measurement program	Lot measurement, time-course measurement (30 s / 60 s intervals)		-
Analysis program	Calibration curve function, concentration conversion		-

Ratings

Item	Absorption method type	Fluorescence method type	Unit
Number of measurement lines	Max.	6	lines
Measurement line spacing (center-to-center spacing)	Min.	3	mm
Measurement line width	0.8 to 1.2		mm
Measurable range	3 × 20		mm
Minimum sensitivity	5	-	mABS ^②
Maximum sensitivity	800	-	mABS ^②
Fluorescence intensity resolution	-	12	bit
Measurement reproducibility ^③	Max.	3 % CV (at +25 °C)	
Operating ambient temperature	+15 to +30		°C
Operating ambient humidity ^④	Bellow 80		% RH
Storage temperature	-20 to +50		°C
Storage humidity ^④	Bellow 80		% RH

② mAbs = milli-absorbance (Absorbance is the logarithm of the ratio of incident to transmitted radiant power through a sample.)

③ When measured using our standard sample. ④ No condensation

Compatible reagent sizes (maximum)

Item	C16723 series	Unit
Reagent housing (W × H × D)	39.5 × 8.0 × 120.4	mm
Dip type reagent strip (W × D)	10 × 120	mm

Lateral flow reader engine C16171 series



The C16171 series lateral flow reader engine consists of an optical head for detecting reagent color development and fluorescence lines and a control board integrated with a motor for driving the optical head. This lateral flow reader engine detects lines with high sensitivity by utilizing software programs for profile data acquisition and waveform analysis. When developing a lateral flow reader, this lateral flow reader engine will serve as the core component specially optimized for reading of lines.


Item	Description / value				
Suffix	-10	-11	-50	-51	-60
Target label	Red-based color lines	Blue-based color lines	Europium	FITC	Europium
Measurement target ①	Color labels		Fluorescent labels		
Measurement method	Absorption method		Fluorescence method		High-sensitivity Fluorescence method
Light source	Green LED	Red LED	UV LED	Blue LED	UV LED

① Two colors cannot be measured simultaneously.

C16171 series operation sequence

When the lateral flow reader engine receives a measurement start command, the optical head starts scanning the reagent color development (fluorescence) area. One measurement completes in about 12 seconds.


• Lateral flow reader engine operation viewed from side



Item	Description / value	Unit
Focal length	6	mm
Self-check run time	Typ. 12	s
Measurement run time	Typ. 12*	

* For C16171-60, 17 s are required.

• Compatible test strip




Item	Description / value	Unit
Detection area (lateral flow direction x width)	20 x 3	mm
Maximum number of detection lines	6	Lines
Detection line width	0.8 to 1.2	mm
Detection line interval	3	

Data acquisition

Reagent color development and fluorescence profile are measured while scan. A measurement is taken a total of 1000 profile data are measured and acquired in one measurement cycle.

• Example of acquired profile (fluorescence line)




Data analysis


Color development and fluorescence lines are automatically analyzed from the acquired profile.

Peak position (Peak-X)
Peak intensity (Peak)
Background (Background)
Area (Area)

• Example of waveform (fluorescence line)



Dimensional outlines (unit: mm)



specifications

Ratings

Item		Description / value		Unit
Supply voltage	Min.	+3.0		V
	Typ.	+3.3 / +5.0		V
	Max.	+5.5		V
Supply current	At power-down	Typ.	0	A
	In standby	Typ.	0.1 ^①	A
	In operation	Max.	1	A
Control input voltage	R x D	Max.	+5.5	V
	PWRDWN	Max.	+5.5	V
Control output voltage	T x D	Max.	+5.5	V
Operating ambient temperature		+15 to +30		°C
Operating ambient relative humidity ^②		Below 80		% RH
Storage temperature		-20 to +50		°C
Storage ambient relative humidity ^②		Below 80		% RH

① C16171-60: 0.2 A

② No condensation

OEM product based on basic model



We offer OEM equipment based on the C16723. A variety of options and accessories can be combined with lateral flow reader engine of core parts, to enable flexible equipment designs.

Basic configuration and functions

• Touch panel



Enables smart and intuitive operation.

• Alarm Function



Alarm sounds when the measurement is finished.

• Memory function



Save the measurement results.

• USB / LAN connection port



Store the measurement data into USB memory sticks.

Optional designs

• GUI design (Software)



GUI customization.

• Sample holder



Customization of holder to fit your reagent kits.

• Panel design



Customization the color of the housing panel and the logo design.

• Wireless communication



Wi-Fi / Bluetooth connectivity can be installed.

Optional accessories

• Code reader



1D or 2D code reading devices.

• External Printer



Prints measurement results.

• Thermal monitor



Monitoring function for ambient temperature.

Product development flow for OEM

01 Inquiry and consultation

After getting a customer inquiry we hold a meeting to ascertain the exact requirements of the customer. Feel free to inquire about anything including device design such as for applications or measurement targets and even device development schedule.




02 Evaluate sample

We provide an evaluation unit for evaluating the sample. You can use this to confirm the reagent compatibility with the evaluation unit, measurement accuracy, and reproducibility. Note that the evaluation unit we provide for making sample evaluations is basically the C16723 series.




03 Evaluate product specifications

We consider the detailed specifications for the device based on results from evaluating the sample. We then define the product specifications according to the reagent measurement results and design a device that maximizes detection performance such as sensitivity and reproducibility.



04 Develop a prototype and test it

We design and fabricate a prototype that meets the specifications in the evaluation. The period for fabricating the prototype device will vary greatly depending on the device specifications and requirements.




05 Make commercial product and mass produce it

We decide on the specifications for the final product (mass-production product). The final product specifications are optimized taking into account the target production volume, material procurement, and sales schedule.




06 Deliver, supply and service

We provide after-follow-up service as needed such as production adjustments and technical support for the finished product.





Accessories and options




Accessories

● Universal holder



- Dedicated AC adapter
Cable length: 1200 mm
- USB cable
- CD-ROM (contains software and user manual)
- Color sample / Fluorescence color sample for checking main unit
Absorption method
C16723-10/-11: A10792
Fluorescence method
C16723-50: A15924
C16723-51: A16441
C16723-60: A16350



Options (sold separately)

● Custom holders

We design and fabricate custom holders that will perfectly match your reagent housing. If our standard Universal holder does not support the size of the reagent housing you are using, our custom-designed holders may prove an ideal solution. So please contact us and provide the reagent housing dimensions.



Custom holder examples

● Dip type holder A10793-01

This holder is designed for our lateral flow readers to allow easy setup of a reagent kit with no housing or a dip type.

Maintenance and inspection

To maintain the reliability and measurement accuracy of the C16723 series lateral flow readers, we recommend making daily checks and doing regular maintenance.

Daily check

You can easily make a daily check by using dedicated a color sample.
Make the check both before and after using the lateral flow reader helps maintain its high reliability.

Recommended check period

- Daily check: Once a day
- Replacement of color sample: Once a year (Fee will be charged)



- Color sample for checking main unit and fluorescence color sample

- Color sample replacement period

The line of the color sample deteriorates due to daily use, so the color sample should be replaced every 12 months. Missing this replacement period may cause discrepancies between the reference value and the actual line density. Do not continue using the color sample if you kept using it after the scheduled replacement period. Purchase a new color sample from us as soon as possible.



Periodic inspections (overhaul)

Please return the unit to us for inspection and maintenance of the entire unit. We will disassemble and inspect individual parts while checking measurement performance and mechanical operation.

Recommended period

- Inspection: Once every 3 years (Fee will be charged)

Frequently Asked Questions (FAQ)

Q. The reagent I'm using does not fit the supplied tray. What should I do?

A. If our Universal holder (see page 9) does not fit your reagent, then our custom-designed holder may be the perfect solution. Please contact us, providing us with the size of the reagent you are using. We can also lend you a demo unit or we can evaluate your sample at our site if you send it to us.

Q. Is it possible to detect color lines other than red and blue?

A. Measurement is possible but the detection sensitivity may be lower than the specifications listed for red and blue colors.

Q. Do you offer any support for inspection and maintenance of the unit?

A. We provide support for maintenance if you send the unit to us. Please contact our sales office. Before sending the unit to us, please clean and disinfect it in view of the needs for bio-safety.

Q. Can the C16723 series make diagnoses in a clinical environment?

A. No, the C16723 series is not a medical device. The C16723 series cannot be used for making diagnoses in a clinical environment. Use it as a tool only for development and management of reagents.

Q. Are there PCs that come with an lateral flow reader?

A. There is no PC that comes with an lateral flow reader. Please prepare one on your own. (We recommend using a Windows® PC with a USB 1.1/2.0 port.)

Q. Is it possible to measure two color lines with one unit?

A. No, multiple color lines cannot be measured with just one unit. A single unit can only measure one color line.

THETA is a registered trademark of Hamamatsu Photonics K.K. (Registered: Japan, EU, U.K., China)

The products listed in this catalog are not finished medical devices and do not have clearance from the FDA or any other medical regulatory body. They are intended as components to be incorporated into OEM systems and it is incumbent on the OEM to obtain a 510(k) clearance for any medical use. These products are not intended for direct sale to dentists, doctors, veterinarians or other healthcare practitioners.

HAMAMATSU PHOTONICS K.K. www.hamamatsu.com

Electron Tube Division

314-5, Shimokanzo, Iwata City, Shizuoka Pref., 438-0193, Japan, Telephone: (81)539/62-5248, Fax: (81)539/62-2205

U.S.A.: HAMAMATSU CORPORATION: 360 Foothill Road, Bridgewater, NJ 08807, U.S.A., Telephone: (1)908-231-0960, Fax: (1)908-231-1218 E-mail: usa@hamamatsu.com

Germany: HAMAMATSU PHOTONICS DEUTSCHLAND GMBH.: Arzbergerstr. 10, 82211 Herrsching am Ammersee, Germany, Telephone: (49)8152-375-0, Fax: (49)8152-265-8 E-mail: info@hamamatsu.de

France: HAMAMATSU PHOTONICS FRANCE S.A.R.L.: 19 Rue du Saule Trapu, Parc du Moulin de Massy, 91882 Massy Cedex, France, Telephone: (33)1 69 53 71 00, Fax: (33)1 69 53 71 10 E-mail: infos@hamamatsu.fr

United Kingdom: HAMAMATSU PHOTONICS UK LIMITED: 2 Howard Court, 10 Twin Road, Welwyn Garden City, Hertfordshire, AL7 1BW, UK, Telephone: (44)1707-294888, Fax: (44)1707-325777 E-mail: info@hamamatsu.co.uk

North Europe: HAMAMATSU PHOTONICS NORDEN AB: Torshamnsgatan 35, 16440 Kista, Sweden, Telephone: (46)8-509-031-00, Fax: (46)8-509-031-01 E-mail: info@hamamatsu.se

Italy: HAMAMATSU PHOTONICS ITALIA S.R.L.: Strada della Moia, 1 int. 6 20044 Arese (Milano), Italy, Telephone: (39)02-93 58 17 33, Fax: (39)02-93 58 17 41 E-mail: info@hamamatsu.it

China: HAMAMATSU PHOTONICS (CHINA) CO., LTD.: 1201, Tower B, Jiaming Center, 27 Dongsanhuan Beilu, Chaoyang District, 100020 Beijing, P.R. China, Telephone: (86)10-6586-6006, Fax: (86)10-6586-2866 E-mail: hpc@hamamatsu.com.cn

Taiwan: HAMAMATSU PHOTONICS TAIWAN CO., LTD.: 13F-1, No.101, Section 2, Gongdao 5th Road, East Dist., Hsinchu City, 300046, Taiwan(R.O.C) Telephone: (886)3-659-0080, Fax: (886)3-659-0081 E-mail: info@hamamatsu.com.tw

TPMZ1031E01
APR. 2024 IP