

NEW



MoxiePlex

Multi-spectral imaging system

C16919-01



The MoxiePlex® Multi-spectral imaging system captures high-resolution images of entire tissue samples stained with multiple fluorescent dyes. These images play a pivotal role in evaluating and visualizing biomarker localization and expression levels, empowering researchers to perform comprehensive analysis.

HAMAMATSU
PHOTON IS OUR BUSINESS


Optimized for fluorescence multiplex analysis imaging system

Acquire high-resolution images with up to 9 channels

The MoxiePlex is designed to speed up multiplex immunofluorescence research by providing up to 9 channels across the entire tissue section.

With a 16 bit output, it delivers high-resolution imaging essential for image analysis. Additionally, the system is equipped with a state-of-the-art camera featuring high quantum efficiency across a 400 nm to 900 nm wavelength range, allowing for a broader selection of fluorescent reagents and enhanced biomarker detection.



Spectral sensitivity



Imaging examples (Human tonsil)



Akoya Biosciences MOTiF PD-1/PD-L1 Panel 7plex kit

Target	Reagent
Nucleus	DAPI
CD8	Opal Polaris 480
PD-L1	Opal 520
FoxP3	Opal 570
PD-1	Opal 620
Cytokeratin	Opal 690
CD68	Opal Polaris 780




Akoya Biosciences MOTiF PD-1/PD-L1 Panel 5plex kit

Target	Reagent
Nucleus	DAPI
PD-L1	Opal 520
FoxP3	Opal 570
Cytokeratin	Opal 690
CD68	Opal Polaris 780




Compatible with BF

H&E



IHC




The MoxiePlex enables seamless superimposition of brightfield and fluorescence images within the Moxie software. It captures high-resolution images of H&E-stained, immunohistochemically (IHC)-stained, and specially stained samples, providing a comprehensive view of tissue morphology for in-depth observation and analysis.


Automated sample recognition and exposure time calculation

The MoxiePlex uses proprietary fluorescence optics technology and algorithms to automate sample recognition, scan range setting, focus position setting, and exposure time calculation. Users can choose between fully automated operation to minimize manual effort or manual adjustments for greater flexibility and control.

Conventional flow




MoxiePlex automatic flow




*3rd party image analysis

• Features of each flow




Automatically recognizes the entire tissue sample and automatically calculates exposure time to acquire fluorescent images

The MoxiePlex enhances immunofluorescence imaging by automatically detecting tissue samples that are challenging to identify in brightfield imaging. It intelligently sets scan profiles and optimizes exposure time for precise fluorescent image acquisition, ensuring a seamless workflow and high accuracy. (Patent pending)



System configuration




Specifications

Product name	MoxiePlex Multi-spectral imaging system	
Product number	C16919-01	
Cassette loader	Up to 60 slides	
Compatible glass slide	75.0 mm to 76.0 mm × 25.0 mm to 26.0 mm (Thickness: 0.9 mm to 1.2 mm)	
Objective lens	20× NA 0.8	
Scanning resolution	20× mode	0.46 μm/pixel
	40× mode	0.23 μm/pixel
Brightfield scan speed	20× mode	Approx. 60 s (15 mm × 15 mm)
	40× mode	Approx. 70 s (15 mm × 15 mm)
Fluorescence throughput *1	20× mode	Approx. 12 min (15 mm × 15 mm, 5 channels)
	40× mode	Approx. 14 min (15 mm × 15 mm, 5 channels)
Fluorescence camera	In-house sCMOS camera	
Z-stack feature	Included	
Image format	OME-TIFF (16 bit, 8 bit), NDPI (8 bit)	
Power supply	AC 100 V to AC 240 V	
Power consumption	Approx. 180 VA	

*1 Fluorescence throughput: The time from loading the glass slide, macro photography, prefocusing 5 points, scanning, image processing, to unloading the glass slide.

Dimensional outlines (Unit: mm)

Main unit: Approx. 95 kg



Filters

Standard compatible reagent	DAPI/FITC/TRITC/Cy5/Cy7/CFP/YFP/mCherry
-----------------------------	---

*Please contact us for details.

- MoxiePlex is a registered trademark of Hamamatsu Photonics K.K. (China, EU, Japan, Switzerland, UK).
 - The product and software package names noted in this brochure are trademarks or registered trademarks of their respective manufacturers.
 - Subject to local technical requirements and regulations, availability of products included in this brochure may vary. Please consult your local sales representative.
 - The product described in this brochure is designed to meet the written specifications, when used strictly in accordance with all instructions.
 - The spectral response specified in this brochure is typical value and not guaranteed.
 - The measurement examples in this brochure are not guaranteed.
 - The products in this brochure are not medical devices.
 - Specifications and external appearance are subject to change without notice.
- © 2025 Hamamatsu Photonics K.K.

HAMAMATSU PHOTONICS K.K. www.hamamatsu.com

Image and Measurement Systems Sales

812, Joko-cho, Chuo-ku, Hamamatsu City, Shizuoka Pref., 431-3196, Japan, Telephone: (81)53-431-0124, Fax: (81)53-433-8031, E-mail: export@sys.hpk.co.jp

U.S.A.: HAMAMATSU CORPORATION: 360 Foothill Road, Bridgewater, NJ 08807, U.S.A., Telephone: (1)908-231-0960, Fax: (1)908-231-1218

Germany: HAMAMATSU PHOTONICS DEUTSCHLAND GMBH: Arzbergerstr. 10, 82211 Herrsching am Ammersee, Germany, Telephone: (49)8152-375-0, Fax: (49)8152-265-8 E-mail: info@hamamatsu.de

France: HAMAMATSU PHOTONICS FRANCE S.A.R.L.: 19 Rue du Saule Trapu, Parc du Moulin de Massy, 91882 Massy Cedex, France, Telephone: (33)1 69 53 71 00, Fax: (33)1 69 53 71 10 E-mail: infos@hamamatsu.fr

United Kingdom: HAMAMATSU PHOTONICS UK LIMITED: 2 Howard Court, 10 Tewin Road, Welwyn Garden City, Hertfordshire, AL7 1BW, UK, Telephone: (44)1707-294888, Fax: (44)1707-325777 E-mail: info@hamamatsu.co.uk

North Europe: HAMAMATSU PHOTONICS NORDEN AB: Torshamnsgatan 35 16440 Kista, Sweden, Telephone: (46)8-509 031 00, Fax: (46)8-509 031 01 E-mail: info@hamamatsu.se

Italy: HAMAMATSU PHOTONICS ITALIA S.R.L.: Strada della Moia, 1 int. 6, 20044 Arese (Milano), Italy, Telephone: (39)02-93 58 17 33, Fax: (39)02-93 58 17 41 E-mail: info@hamamatsu.it

China: HAMAMATSU PHOTONICS (CHINA) CO., LTD.: 1201 Tower B, Jiaming Center, 27 Dongsanhuan Bellu, Chaoyang District, 100020 Beijing, P.R. China, Telephone: (86)10-6586-6006, Fax: (86)10-6586-2866 E-mail: hpc@hamamatsu.com.cn

Taiwan: HAMAMATSU PHOTONICS TAIWAN CO., LTD.: 13F-1, No.101, Section 2, Gongdao 5th Road, East Dist., Hsinchu City, 300046, Taiwan(R.O.C.), Telephone: (886)3-659-0080, Fax: (886)3-659-0081 E-mail: info@hamamatsu.com.tw

Korea: HAMAMATSU PHOTONICS KOREA CO., LTD.: A-912, 167, Songpa-daero, Seoul, 05855, Republic of Korea, Telephone: (82)2-2054-8202, Fax: (82)2-2054-8207 E-mail: sales@hpk.co.kr

Cat. No. SBIS0152E02
JUN/2025 HPK
Created in Japan