

RUO model

NanoZoomer[®] S540

Digital slide scanner

C17400-01

High throughput scanner with rack compatibility and continuous loading



High-speed scanning

Approx. 30 s/slide
(20×/40× mode)

High throughput

82 slides/h
(20×/40× mode)

Slide rack compatibility

540 slides
Rack for 30 slides
360 slides
Rack for 20 slides

Continuous loading

Load racks
without interrupting
scanning

Intuitive touchscreen interface

Easy to use
interface for
efficient operation

A stable system for operational efficiency

The NanoZoomer S540 has featuring compatibility with staining racks, fast slide scanning, exceptional image quality, and an intuitive touchscreen interface, this high-capacity scanner to handle your needs.

Up to
540
Slides*


82
Slides/h
(20×/40× mode)

*Maximum capacity: 540 slides with a rack for 30 slides and 360 slides with a rack for 20 slides.



Rack compatibility

The NanoZoomer S540 can handle both racks for 30 slides and 20 slides simultaneously. Before installation, users can select a rack configuration with either a 30-slide or 20-slide capacity, designed to fit into 18 slots.



Selectable scanning modes

Fully-Automatic scanning

All scanning processes can be fully automated. Scanning begins immediately upon loading the rack, and the scanner's touchscreen interface provides a clear visual indication of scan progress.

Semi-Automatic scanning

Set-up scanning conditions such as the scan area, focus point placement, Z-stacking, and resolution and to assign bespoke profiles to each slide or rack.

Continuous loading

Insert and exchange racks without interrupting scanning. Rack priority can be set either directly in the touchscreen interface or using the operation software (NZAcquire or NZConnect Scan*).

*Please refer to the datasheet of NZConnect Scan to see the full list of features.



Touchscreen interface

The touchscreen allows users to operate and monitor key scanning functions without accessing the operation software, including:

- Start/Stop scanning and remove racks
- Set scan priority for racks
- Display scan status



Start/Stop scanning and remove racks



Set scan priority for racks



Display scan status

AI tissue detection

AI tissue detection enables fast and targeted tissue identification, ensuring that only tissue areas are scanned accurately.

Profile functions

Scan profile

Save scan conditions such as scan area and sample size threshold for tissue recognition.


Output profile

Assign output file conditions, such as save location and file type parameters.

Quality check


QC (Quality check) mode is available to allow users to check image quality before accepting the whole slide imaging.

System configuration




Dimensional outlines (Unit: mm)

- Scanner: Approx. 125 kg
- Dedicated cart: Approx. 65 kg



LASER SAFETY

The NanoZoomer S540 is a Class 2 laser product. Hamamatsu Photonics classifies laser diodes and provides appropriate safety measures and labels according to the classification required for manufacturers by IEC 60825-1. When using this product, follow all safety measures as specified by the IEC.



Caution Label

Specifications

Product name		NanoZoomer S540 Digital slide scanner
Product number		C17400-01
Scanning speed*1 (20× mode/40× mode)		Approx. 30 s/slide
Throughput*1 (20× mode/40× mode)		More than 82 slides/h
Objective lens		20× NA 0.75 User can select 20× or 40× mode at start of scanning
Compatible glass slides		25.0 mm to 26.0 mm × 75.0 mm to 76.0 mm (Thickness 0.9 mm to 1.2 mm)
Slide loader	Rack for 20 slides	Max. 360 slides (20 slides × 18 racks)
	Rack for 30 slides*2	Max. 540 slides (30 slides × 18 racks)
Scanning resolution	20× mode	Approx. 0.46 μm/pixel
	40× mode	Approx. 0.23 μm/pixel
Focusing method		Pre-Focus map
Z-stack feature		Included
Image compression		JPEG compression
Readable barcodes	1D Barcodes	Code 39, Code 128, Interleaved 2 of 5, Codabar, EAN-8 and UPC-E
	2D Barcodes	DataMatrix (ECC200) QR Code® (QR Code Model-2)
Power supply		AC 100 V to AC 240 V
Power consumption (Scanner only)		Approx. 200 VA

*1 When scanning, it scans an area of 15 mm × 15 mm square with 5 focus points.
*2 This configuration is provided with the rack for 20 slides.

- NanoZoomer and NDP are registered trademarks of Hamamatsu Photonics K.K. (EU, Japan, UK, USA).
 - QR Code is a registered trademark of DENSO WAVE INCORPORATED.
 - The product and software package names noted in this brochure are trademarks or registered trademarks of their respective manufacturers.
 - Subject to local technical requirements and regulations, availability of products included in this brochure may vary. Please consult your local sales representative.
 - The product described in this brochure is designed to meet the written specifications, when used strictly in accordance with all instructions.
 - The products in this brochure are not medical devices.
 - Specifications and external appearance are subject to change without notice.
- ©2025 Hamamatsu Photonics K.K.

HAMAMATSU PHOTONICS K.K. www.hamamatsu.com

Image and Measurement Systems Sales
812, Joko-cho, Chuo-ku, Hamamatsu City, Shizuoka Pref., 431-3196, Japan, Telephone: (81)53-431-0124, Fax: (81)53-433-8031, E-mail: export@sys.hpk.co.jp

U.S.A.: HAMAMATSU CORPORATION: 360 Foothill Road, Bridgewater, NJ 08807, U.S.A., Telephone: (1)908-231-0960, Fax: (1)908-231-1218
Germany: HAMAMATSU PHOTONICS DEUTSCHLAND GMBH: Arzbergerstr. 10, 82211 Herrsching am Ammersee, Germany, Telephone: (49)8152-375-0, Fax: (49)8152-265-8 E-mail: info@hamamatsu.de
France: HAMAMATSU PHOTONICS FRANCE S.A.R.L.: 19 Rue du Saule Trapu, Parc du Moulin de Massy, 91882 Massy Cedex, France, Telephone: (33)1 69 53 71 00, Fax: (33)1 69 53 71 10 E-mail: infos@hamamatsu.fr
United Kingdom: HAMAMATSU PHOTONICS UK LIMITED: 2 Howard Court, 10 Tewin Road, Welwyn Garden City, Hertfordshire, AL7 1BW, UK, Telephone: (44)1707-294888, Fax: (44)1707-325777 E-mail: info@hamamatsu.co.uk
North Europe: HAMAMATSU PHOTONICS NORDEN AB: Torshamnsgatan 35 16440 Kista, Sweden, Telephone: (46)8-509 031 00, Fax: (46)8-509 031 01 E-mail: info@hamamatsu.se
Italy: HAMAMATSU PHOTONICS ITALIA S.R.L.: Strada della Moia, 1 int. 6, 20044 Arese (Milano), Italy, Telephone: (39)02-93 58 17 33, Fax: (39)02-93 58 17 41 E-mail: info@hamamatsu.it
China: HAMAMATSU PHOTONICS (CHINA) CO., LTD.: 1201 Tower B, Jiaming Center, 27 Dongsanhuan Bellu, Chaoyang District, 100020 Beijing, P.R. China, Telephone: (86)10-6586-6006, Fax: (86)10-6586-2866 E-mail: hpc@hamamatsu.com.cn
Taiwan: HAMAMATSU PHOTONICS TAIWAN CO., LTD.: 13F-1, No.101, Section 2, Gongdao 5th Road, East Dist., Hsinchu City, 300046, Taiwan(R.O.C.), Telephone: (886)3-659-0080, Fax: (886)3-659-0081 E-mail: info@hamamatsu.com.tw
Korea: HAMAMATSU PHOTONICS KOREA CO., LTD.: A-912, 167, Songpa-daero, Seoul, 05855, Republic of Korea, Telephone: (82)2-2054-8202, Fax: (82)2-2054-8207 E-mail: sales@hpk.co.kr