


CAMERA LINE UP CATALOG



APPLICATION



Application & Case study


https://camera.hamamatsu.com/all/en/application_and_case_study.html

We have a diverse lineup of cameras that support a wide range of wavelengths from X-rays to the near-infrared and support a variety of applications.

Life science

Super resolution microscopy


Cells are observed with higher spatial resolution than diffraction limit by the super resolution microscopy.



Camera: ORCA®-Quest
Super resolution imaging system: VT-iSIM
Data courtesy of Steven Coleman (Visitech International Ltd.)

Bioluminescence measurements


Ligand-stimulated binding of NanoLuc-Arrestin to GPCR-mVenus and its intracellular uptake are observed by simultaneous two-wavelength luminescence imaging.



Data courtesy of Dr. Masataka Yanagawa (Tohoku university)

Light sheet microscope


The zebrafish larvae brain function during its natural behavior is observed with light sheet fluorescence microscope.



Data courtesy of Dr. Drew Robson (Max Planck Institute for Biological Cybernetics)

Observation of cultured cells

Cells cultured in one well of a microplate are observed by high-resolution imaging with fluorescent images.




* Displayed with pseudo color by image processing.

Quantum Technology

Quantum computing (Neutral atom, Ion)


Position and quantum state of Rb atoms, trapped and arranged one by one in vacuum, are diagnosed via fluorescence.



Camera: ORCA-Quest
Data courtesy of Prof. Takashi Yamamoto and Associate Prof. Toshiki Kobayashi (Osaka University)

Quantum optics

qCMOS® camera is used for absorption imaging with quantum light source to compare between 1 photon event and 2 photon event images.




Camera: ORCA-Quest
Data courtesy of Prof. Miles Padgett (University of Glasgow)

Astronomy

Lucky imaging

Wide field of view and low-noise imaging is used to obtain a clear image of the stars by integrating, from among many acquired images, that are less affected by atmospheric turbulence.




* Displayed with pseudo color by image processing.

Synchrotron imaging

X-ray phase contrast CT image of mouse embryo

The mouse embryo is observed using the synchrotron X-ray.




ORCA-Quest combined with High resolution X-ray imaging system (M11427)
Data courtesy of SPring-8 BL20B2 beamline by Dr. Masato Hoshino, Senior Scientist in Japan Synchrotron Radiation Research Institute (JASRI)

Semiconductor inspection

Transmission observation of Si wafer

The pattern formed on the Si wafer is transmissively observed from the backside.




Semiconductor device observation

The pattern under the Si layer is observed by infrared imaging.



Structure observation of semiconductor devices

The interior structure of a semiconductor device is analyzed at the nano-level by high-resolution imaging using an electron microscope.




20 nm

Food inspection

Foreign object detection


Small stones mixed in coffee beans that are difficult to see with visible light is detected by the infrared imaging.



Analysis / Spectroscopy

Material identification

Infrared imaging identifies materials that are difficult to distinguish in visible light, such as PVC, acrylic, PET, and PS.




























* Displayed with pseudo color by image processing.

CAMERA LINE UP



For detailed information

https://www.hamamatsu.com/all/en/product/cameras.html


| Wavelength range | Visible to near-infrared | | | | | | | | | Near-infrared | | | | | Wavelength range | |
|--|---|---|---|--|--|---|---|---|---|---|---|---|---|---|--|--|
| Name | ORCA-Quest 2 qCMOS camera | ORCA-Fire Digital CMOS camera | ORCA-Fusion BT Digital CMOS camera | ORCA-Fusion Digital CMOS camera | ORCA-Flash4.0 V3 Digital CMOS camera | ORCA-Flash4.0 LT3 Digital CMOS camera | ORCA-Halo sCMOS camera | ORCA-spark Digital CMOS camera | TDI camera | InGaAs camera | | | InGaAs line scan camera | Name | | |
| Type | C15550-22UP | C16240-20UP/-30UP | C15440-20UP | C14440-20UP | C13440-20CU | C11440-42U40 | C17440-20U | C11440-36U | C10000-801 | C16741-40U | C14041-10U-02 | C12741-03-02 | C12741-11 | C15333-10E04-02 | Type | |
| Appearance |  |  |  |  |  |  |  |  |  |  |  |  |  |  | Appearance | |
| Image sensor type | Area sensor | | | | | | | | | TDI sensor | Area sensor | | | Line sensor | Image sensor type | |
| Sensitivity wavelength range (nm) (Spectral response: See P5) | 250 to 1000 | | | 350 to 1000 | | | | | | 200 to 1000 | 400 to 1700 | 950 to 1700 | | 900 to 1550 | 950 to 1700 | Sensitivity wavelength range (nm) (Spectral response: See P5) |
| |  |  |  |  |  | |  |  |  | *3 |  | |  |  | | |
| Effective number of pixels (H×V) | 4096 × 2304 | 4432 × 2368 | 2304 × 2304 | 2304 × 2304 | 2048 × 2048 | 2048 × 2048 | 3000×3000 | 1920 × 1200 | 2048 × 128 | 1280 × 1024 | 320 × 256 | 640 × 512 | 640 × 512 | 1024 × 1 | Effective number of pixels (H×V) | |
| Pixel size ((H) μm × (V) μm) | 4.6 × 4.6 | 4.6 × 4.6 | 6.5 × 6.5 | 6.5 × 6.5 | 6.5 × 6.5 | 6.5 × 6.5 | 3.76×3.76 | 5.86 × 5.86 | 12 × 12 | 5 × 5 | 20 × 20 | | 12.5 × 12.5 | | Pixel size ((H) μm × (V) μm) | |
| Effective area ((H) mm × (V) mm) | 18.841 × 10.598 | 20.387 × 10.892 | 14.976 × 14.976 | 14.976 × 14.976 | 13.312 × 13.312 | 13.312 × 13.312 | 11.280×11.280 | 11.25 × 7.03 | 24.58 × 1.536 | 6.40 × 5.12 | 6.4 × 5.12 | 12.8 × 10.24 | 12.8 × 10.24 | 12.8 × 0.0125 | Effective area ((H) mm × (V) mm) | |
| Full well capacity (electrons) typ.*1 | 7000 | 20 000 | 15 000 | 15 000 | 30 000 | 30 000 | 49 100 | 33 000 | 80 000 | - | - | | 300 000 | - | Full well capacity (electrons) typ.*1 | |
| Dynamic range typ.*1 | 23 000:1 | 20 000:1 | 21 400:1 | 21 400:1 | 37 000:1 | 33 000:1 | 12 000:1 | 5000:1 | 1600:1 | - | - | | - | - | Dynamic range typ.*1 | |
| Cooling method | Forced-air cooled Water cooled | Forced-air cooled (-20UP) Water cooled (-30UP) | Forced-air cooled Water cooled | Forced-air cooled Water cooled | Forced-air cooled Water cooled | Forces-air cooled | Forced-air cooled Water cooled | - | - | Forced-air cooled Natural-air cooled | Forced-air cooled | | Forced-air cooled Water cooled | - | Cooling method | |
| Cooling temperature (°C)*1 | -20 -35 (max cooling) | +20 | -8 -15 (max cooling) | -5 -15 (max cooling) | -10 -30 (max cooling) | +10 | +10 | - | - | +15 (Forced-air cooled) | +10 | | -70 (Water cooled) -60 (Forced-air cooled) | - | Cooling temperature (°C)*1 | |
| Readout speed (frame/s) (Full resolution)*1 | 120 | 115 | 89.1 | 89.1 | 100 | 40 | 24.3 | 64.9 | 50 kHz (Line rate) | 71.53 | 216.6 | 59.774 | 7.2 | 40 kHz (Line rate) | Readout speed (frame/s) (Full resolution)*1 | |
| Readout noise (electrons) rms typ.*1 | 0.30 | 1.0 | 0.7 | 0.7 | 1.4 | 1.5 | 1.3 | 6.6 | 50 | - | - | | 500 | - | Readout noise (electrons) rms typ.*1 | |
| Dark current (electrons/pixel/s) typ.*1 | 0.016 0.006 (max cooling) | 0.6 | 1.0 0.7 (max cooling) | 0.5 0.2 (max cooling) | 0.6 0.006 (max cooling) | 0.6 | 0.03 | - | - | - | - | | 130 (Water cooled) 300 (Forced-air cooled) | - | Dark current (electrons/pixel/s) typ.*1 | |
| Interface | CoaXPress (Quad CXP-6) USB 3.1 Gen1 | CoaXPress (Quad CXP-6) USB 3.1 Gen1 | CoaXPress (Dual CXP-6) USB 3.0 *2 | CoaXPress (Dual CXP-6) USB 3.0 *2 | Camera Link USB 3.0 *2 | USB 3.1 Gen 1 | USB 3.1 Gen1 | USB 3.0 *2 | Camera Link | USB 3.1 Gen 1 | USB 3.0 *2 | USB 3.0 *2/EIA | Camera Link | Gigabit Ethernet | Interface | |
| Applications | Life science imaging Quantum technology Astronomy Semiconductor inspection Synchrotron imaging Electron microscope | Life science imaging Synchrotron imaging Electron microscope | Life science imaging Synchrotron imaging Electron microscope | Life science imaging Semiconductor inspection Synchrotron imaging Electron microscope | Life science imaging Semiconductor inspection Synchrotron imaging Electron microscope | Life science imaging Semiconductor inspection | Life science imaging Synchrotron imaging | Life science imaging Synchrotron imaging | Life science imaging Semiconductor inspection | Semiconductor inspection Food inspection Analysis/spectroscopy | Semiconductor inspection Food inspection Analysis/spectroscopy | | Life science imaging Semiconductor inspection | Semiconductor inspection Food inspection | Applications | |

*1 Depends on the mode and conditions. For details, please refer to each product catalog. *2 Equivalation to USB 3.1 Gen1 *3 For detailed information, please refer to product catalog.

| Camera type | Board type camera for OEM | | | | | For X-ray |
|---|------------------------------------|---------------------------------|-------------|-------------|------------------------|--|
| Name | Scientific CMOS board level camera | Digital CMOS board level camera | | | TDI board level camera | X-ray sCMOS camera |
| Type | C11440-52U30 | C13949-50U | C13770-50U | C13752-50U | C10000-A01 | C12849-111U |
| Appearance | | | | | | |
| Image sensor type | Area sensor | Area sensor | | | TDI sensor | Area sensor |
| Sensitivity wavelength range (nm) (Spectral response: See P5) | 350 to 1000 | 350 to 1000 | | | 200 to 1000 | 25 kV to 90 kV (Recommended X-ray tube voltage range) |
| Effective number of pixels (H × V) | 2048 × 2048 | 4096 × 3008 | 2464 × 2056 | 2048 × 1544 | 2048 × 128 | 2048 × 2048 |
| Pixel size ((H) μm × (V) μm) | 6.5 × 6.5 | | 3.45 × 3.45 | | 12 × 12 | 6.5 × 6.5 |
| Effective area ((H) mm × (V) mm) | 13.312 × 13.312 | 14.13 × 10.37 | 8.50 × 7.09 | 7.06 × 5.32 | 24.53 × 1.536 | 13.312 × 13.312 |
| Full well capacity (electrons) typ.*1 | 30 000 | 10 500 | | | 80 000 | 30 000 |
| Dynamic range typ.*1 | 18 000:1 | 4565:1 | | | 1600:1 | 18 000:1 |
| Readout speed (frame/s) (Full resolution)*1 | 30 | 15 | 40 | 65 | 50 kHz (Line rate) | 30 |
| Readout noise (electrons) rms typ.*1 | 2.3 | 2.3 | | | 50 | 2.3 |
| Interface | USB 3.0 *2 | USB 3.0 *2 | | | Camera Link | USB 3.0 *2 |
| Applications | Contact us | Contact us | Contact us | Contact us | Contact us | Synchrotron imaging |

*1 Depends on the mode and conditions. For details, please refer to each product catalog. *2 Equivalation to USB 3.1 Gen1


Spectral response



OEM CAMERA

We design and manufacture OEM cameras specific to each customer. We provide various types of cameras with options such as shape, sensor, interface, cooling method, software, etc. to meet customers' requests. The measurement wavelength range covers not only the visible range but spans widely from X-ray to infrared.

Cost reduction with minimum required functions




Shorten delivery time with simulation technology

We can perform imaging simulations that match the characteristics of various cameras (wavelength, sensitivity, speed, etc.). By using this technology, we can shorten the process of repeating design and trial production, and provide cameras that meet your purpose efficiently and in a short time.

Flow from inquiry to delivery


Without simulation



With simulation (in our company)




This is a simulation example using a cytoskeleton sample. The amount of light per pixel is set to the same value for simulation. Simulations can be performed by flexibly changing the acquisition conditions such as exposure time, and the results can be viewed not only as still images but also as movies.



You can try the simulation on our website below.
Camera Simulation Engine URL: https://camera.hamamatsu.com/all/en/learn/camera_simulation_engine.html

SOFTWARE

We provide a common camera library "DCAM-API," Hamamatsu Photonics software that can maximize the characteristics of your camera, and a tool "DCAM-SDK," that allows you to build your own control software. Through DCAM-API, even if the camera or interface is changed, the software modification/change can be minimized.



Third-party software

Our cameras can be controlled by combining our cameras and peripherals with software from each microscope manufacturer, bioimaging software, or the following software.

Plugins that are compatible with third-party software

| Software | Manufacturer | OS |
|----------|----------------------|---------|
| LabView | National Instruments | Windows |
| MATLAB® | The MathWorks | Windows |
| µManager | Open source | Windows |
| EPICS | Open source | Linux |

*For details on external software, please contact the manufacturer.

Please download plugins from the URL below.
<https://dcam-api.com/third-party-plugins/>

For details, please refer to the following link. <https://www.hamamatsu.com/all/en/product/cameras/software/third-party-software.html>



RELATED PRODUCTS

Imaging optical system

We also have a lineup of Imaging optical systems to expand usability of our cameras, such as multi wavelength imaging and High resolution X-ray imaging system.



Image splitting optics W-VIEW GEMINI A12801

Product details page URL:
<https://www.hamamatsu.com/all/en/product/optical-components/image-splitting-optics.html>



High resolution X-ray imaging system M11427

Product details page URL:
<https://www.hamamatsu.com/all/en/product/cameras/high-resolution-x-ray-imaging-system.html>

X-ray line scan camera/X-ray TDI camera

We have a lineup of X-ray non-destructive inspection cameras that can be used in-line. Since it is possible to inspect the inside of substances that cannot be seen with visible light or infrared light, these cameras are suitable for foreign matter inspection of foods and pharmaceuticals, defect inspection of printed circuit board, etc.




X-ray line scan camera C14300 series

Product details page URL:
<https://www.hamamatsu.com/all/en/product/cameras/x-ray-line-scan-cameras/C14300series.html>



X-ray TDI camera C12300 series

Product details page URL:
<https://www.hamamatsu.com/all/en/product/cameras/x-ray-tdi-cameras.html>



- DCAM-API is a registered trademark of Hamamatsu Photonics K.K. (EU, Japan, UK, USA).
 - ORCA and qCMOS are registered trademarks of Hamamatsu Photonics K.K. (China, EU, Japan, UK, USA).
 - Windows is a registered trademark of Microsoft Corporation in the USA (and other countries).
 - MATLAB is a registered trademark of the MathWorks, Inc.
 - The product and software package names noted in this brochure are trademarks or registered trademarks of their respective manufacturers.
 - Subject to local technical requirements and regulations, availability of products included in this brochure may vary. Please consult your local sales representative.
 - The products described in this brochure are designed to meet the written specifications, when used strictly in accordance with all instructions.
 - The spectral response specified in this brochure is typical value and not guaranteed.
 - Specifications and external appearance are subject to change without notice.
- © 2025 Hamamatsu Photonics K.K.

HAMAMATSU PHOTONICS K.K. www.hamamatsu.com

Image and Measurement Systems Sales

812, Joko-cho, Chuo-ku, Hamamatsu City, Shizuoka Pref., 431-3196, Japan, Telephone: (81)53-431-0124, Fax: (81)53-433-8031, E-mail: export@sys.hpk.co.jp

U.S.A.: HAMAMATSU CORPORATION: 360 Foothill Road, Bridgewater, NJ 08807, U.S.A., Telephone: (1)908-231-0960, Fax: (1)908-231-1218

Germany: HAMAMATSU PHOTONICS DEUTSCHLAND GMBH: Arzbergerstr. 10, 82211 Herrsching am Ammersee, Germany, Telephone: (49)8152-375-0, Fax: (49)8152-265-8 E-mail: info@hamamatsu.de

France: HAMAMATSU PHOTONICS FRANCE S.A.R.L.: 19 Rue du Saule Trapu, Parc du Moulin de Massy, 91882 Massy Cedex, France, Telephone: (33)1 69 53 71 00, Fax: (33)1 69 53 71 10 E-mail: infos@hamamatsu.fr

United Kingdom: HAMAMATSU PHOTONICS UK LIMITED: 2 Howard Court, 10 Tewin Road, Welwyn Garden City, Hertfordshire, AL7 1BW, UK, Telephone: (44)1707-294888, Fax: (44)1707-325777 E-mail: info@hamamatsu.co.uk

North Europe: HAMAMATSU PHOTONICS NORDEN AB: Torshamnsgatan 35 16440 Kista, Sweden, Telephone: (46)8-509 031 00, Fax: (46)8-509 031 01 E-mail: info@hamamatsu.se

Italy: HAMAMATSU PHOTONICS ITALIA S.R.L.: Strada della Moia, 1 int. 6, 20044 Arese (Milano), Italy, Telephone: (39)02-93 58 17 33, Fax: (39)02-93 58 17 41 E-mail: info@hamamatsu.it

China: HAMAMATSU PHOTONICS (CHINA) CO., LTD.: 1201 Tower B, Jiaoming Center, 27 Dongsanhuan Bellu, Chaoyang District, 100020 Beijing, P.R. China, Telephone: (86)10-6586-6006, Fax: (86)10-6586-2866 E-mail: hpc@hamamatsu.com.cn

Taiwan: HAMAMATSU PHOTONICS TAIWAN CO., LTD.: 13F-1, No.101, Section 2, Gongdao 5th Road, East Dist., Hsinchu City, 300046, Taiwan(R.O.C.), Telephone: (886)3-659-0080, Fax: (886)3-659-0081 E-mail: info@hamamatsu.com.tw

Korea: HAMAMATSU PHOTONICS KOREA CO., LTD.: A-912, 167, Songpa-daero, Seoul, 05855, Republic of Korea, Telephone: (82)2-2054-8202, Fax: (82)2-2054-8207 E-mail: sales@hpk.co.kr

Cat. No. SCAS0142E06
JUN/2025 HPK
Created in Japan