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Poropare™

取扱説明書 (DESI 用) / User Manual (for DESI) Ver. 1.1

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Table of contents (English edition)

1 Notes on using this product	9
2 Safety precautions	10
3 Handling precautions	10
4 Product description	11
5 Handling method	11
5.1 Unpacking and storage	11
5.2 Taking the transfer plate out of the case	12
5.3 Transferring the sample	12
5.4 Attaching to the adapter	12
6 Disposal	12
7 When you are in need of help	12

1 Notes on using this product

- Before using this product be sure to read the latest version of the user manual (this manual) and follow the safety precautions to ensure correct use.
 - Keep this product out of the reach of persons not familiar with this product.
 - We make constant efforts to improve the product quality and reliability but this does not guarantee the complete integrity of the product.
 - This product was developed and manufactured to assist in sample preparation for mass spectrometry and its imaging by utilizing the desorption electrospray ionization (DESI)
 - To obtain detailed information on how to install this product in your mass spectrometer and how to use it after installation, please contact the dealer where you purchased this product or consult the mass spectrometer manufacturer.
 - We cannot be held liable for any damage caused by natural disasters such as earthquakes, lightning, windstorms, and floods, or acts of third parties, other accidents, intentional or negligent misuse by the customer, or usage under other unusual conditions.
 - We assume no responsibility for any incidental damage (failure of peripheral equipment including mass spectrometers, loss of business profits, or business interruptions) arising from the use or inability to use this product.
 - Maintenance of peripheral equipment including mass spectrometers and replacement of consumable parts associated with the use of this product shall be performed by the customer on their own responsibility. Please note that mass spectrometry imaging measurements generally cause wear in mass spectrometers more rapidly than in ordinary spectrometry.
 - Measurement results obtained using this product will vary depending on the method of use and the operating environment. This product is intended to assist in sample preparation and does not guarantee that the desired measurement results will be obtained.
 - Although every effort is made to ensure product quality, in the case of product damage caused during shipping or discoloration or deterioration (which affects measurement) found immediately after unpacking, we will replace the product with a new one, provided that all enclosed transfer plates are unused. In such a case, please contact the dealer where you purchased this product and then return the product along with the package described the lot number. If no complaint is reported within 3 months from first receiving this product, the product will be deemed to have been accepted by the user as a good product.
-
- ✧ Reproduction or copying of this user manual in whole or in part is prohibited without our express permission.
 - ✧ The contents of this user manual are subject to change without prior notice.
 - ✧ The contents of this user manual are believed to be reliable. However, if you find any part of it unclear, then please contact us for assistance.

2 Safety precautions

Failure to observe the following precautions may result in injuries, equipment malfunctions or other unforeseen accidents.


- Do not use this product for samples that are generally considered hazardous.
Failure to do so may cause fire, burns, injuries, or equipment malfunctions. This product is not intended for use with samples that are explosive, ignitable, flammable, or otherwise hazardous.
- Do not apply excessive force to this product. Also, do not use this product under conditions where a localized force is applied (for example, due to using a pointed hard sample).
This product is made of glass so there is a risk of breakage.
- If this product is damaged, do not touch it with bare hands.
The broken surface is sharp and might cause injuries from cuts.
- Do not install this product in the equipment if damaged.
This product might otherwise drop off into the equipment or cause damage to the equipment.
- Be sure to use the appropriate adapter for this product and your system.
This product might otherwise drop off into the equipment or cause damage to the equipment.
- Securely install this product in place.
This product might otherwise drop off into the equipment or cause damage to the equipment.

3 Handling precautions


Failure to observe the following precautions may prevent this product from delivering full performance or might adversely affect the equipment.

- Avoid handling this product with bare hands and do not touch the active surface.
Unintended signals might be detected due to foreign matter or dirt that adheres to the product. When using gloves, use powder-free gloves.
- Do not reuse this product.
This product is disposable and cannot be reused. Reuse may generate unintended signals.
- Do not transfer highly adhesive samples or fragile, crumbly samples.
Excessive sample amounts adhering to the active surface of this product might cause unexpected problems such as sample scattering during measurement.
- Store this product in a low-humidity environment and do not expose it to normal atmospheric conditions for long periods of time.
The active surface might lose its function or cause unintended signals to appear.

4 Product description

Appearance

Active area detailed view


Materials used	Soda glass
Approximate product dimensions	76 mm on long side, 26 mm on short side, 1.05 mm (max.) in plate thickness
Maximum effective area	71 mm on long side, 21 mm on short side
Usage environment	Room temperature range (+10 °C to +30 °C)

5 Handling method

5.1 Unpacking and storage

Packaged state


Case enclosed in package
(Effective surface facing upward)

When unpacking, be careful not to drop the case inside the package. This packaging is moisture-proofed to keep the product in a suitable environment during shipping. Please use the product as soon as possible after delivery. Be aware that once the package is opened, the moisture-proof effect will be lost, so open the package just before use. Please note the direction of the effective surface since the transfer plate is enclosed in the case as shown above. Since the case has no moisture-proofing, do not leave it in a room environment after taking it out of the package, and be sure to perform the sample transfer as soon as possible. When storing this product after unpacking, we recommend storing it in a hermetically sealed environment filled with a low-moisture inert gas such as in a dehumidified nitrogen desiccator. However, this is not a guarantee of product performance after unpacking and storing the product.

5.2 Taking the transfer plate out of the case

Avoid handling the transfer plate including the edges with bare hands and take care not to touch the active surface. When taking a transfer plate out of the case, first check which side is the active surface and then, while holding the transfer plate by the edges with fingers, take it out from the slits holding it. At this point, be careful to avoid applying excessive force to the transfer plate.

Wear powder-free gloves and take out the transfer plate by pulling straight up from the slits while being careful not to touch the active surface.

5.3 Transferring the sample

When transferring a sample, always comply with the instructions in "2. Safety precautions" and "3. Handling precautions." To prevent the transfer plate from bending or warping, place it on a flat hard surface so that the entire surface of the transfer plate is in uniform contact with the flat hard surface. However, please be careful to avoid cracking or damaging the transfer plate during the sample transfer. After transfer, check that there are no dust particles or foreign matter on the surface, and remove them if found.

5.4 Attaching to the adapter

Before installing the transfer plate in the equipment, be sure to remove any foreign matter and dust particles by using an air duster or similar device and be sure to check the transfer plate for any damage.

Use the appropriate adapter for this product and your instrument and be installed to the adapter without loosening. When the product does not fit into the adapter, do not force it, but contact the manufacturer. Set the measurement conditions according to the instructions in the mass spectrometer user manual.

6 Disposal

When disposing of this product, take appropriate measures in compliance with applicable regulations regarding waste disposal and correctly dispose of it yourself or entrust proper disposal to a licensed industrial waste disposal company. In any case, be sure to comply with the regulations in your country or state to ensure correct disposal.

7 When you are in need of help

To find more information on Poropare™ transfer plates, please access our website (<https://www.hamamatsu.com/all/en/product/optical-components/transfer-plate.html>)

There you can see the related information and refer to it as needed.

To find more detailed information on mass spectrometers, please contact the equipment manufacturer.

Reproduction or copying of this user manual in whole or in part is prohibited without our express permission.

The contents of this user manual are believed to be reliable. However, if you find any part of it unclear, then please contact us for assistance.

The contents of this user manual are subject to change without prior notice.

HAMAMATSU PHOTONICS K.K. www.hamamatsu.com

Electron Tube Division
314-5, Shimokanzo, Iwata City, Shizuoka Pref., 438-0193, Japan, Telephone: (81)539/62-5248, Fax: (81)539/62-2205

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