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H17362 Series

取扱説明書 / Instruction manual

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Table of contents

1. Safety precautions	20
1-1. Alert symbols and signal words	20
1-2. Safety precautions	20
1-3. Disposal	21
1-4. Registered trademarks	21
1-5. FCC statement.....	22
2. Overview	23
3. Part names	24
4. Package contents	25
5. Setup	25
5-1. Installing the device driver	25
5-2. Setting up the sample software	26
6. How to use the sample software	27
6-1. Sample software screen	27
6-2. Operating procedures	28
7. Specifications.....	29
7-1. General ratings	29
7-2. Specifications.....	29
7-3. Other Specifications.....	30
7-4. Photodetector characteristics	31
7-5. Dimensional outlines.....	33
7-6. Block diagram	34
7-7. DLL file instructions.....	35
8. Warranty	38
9. Repair	38

1. Safety precautions

1-1. Alert symbols and signal words

To help understand the safety instructions in this manual, the following symbols and signal words are used. Make sure that you understand their contents and comply with the instructions.

WARNING

"WARNING" indicates a potentially hazardous situation which, if not avoided, could result in death or serious injury.

CAUTION

"CAUTION" indicates a potentially hazardous situation which, if not avoided, could result in minor or moderate injury or damage to the equipment.

This symbol shows a caution or warning you must observe. Read the instructions carefully to ensure correct and safe product use.

This symbol indicates a prohibited action. Read the instructions carefully and never attempt the prohibited action

This symbol indicates a mandatory action. Read the instructions carefully and always follow them.

1-2. Safety precautions

WARNING

Do not remove any cover of this unit and do not disassemble this unit. Doing so may damage this unit or cause electrical shock since some internal parts generate a high voltage.

Do not handle this unit with wet hands. Doing so may cause electrical shock. This unit is not designed to be waterproof, so do not immerse it in water.

CAUTION

Some parts of the photodetector in this unit are made of glass, so excessive shock or force may cause cracks. Do not apply excessive shock or force to this unit.

This unit is not designed for use in a vacuum or in a corrosive gas atmosphere. Do not use this unit in such environments.

Do not use this unit in humid environments where condensation occurs. Condensation may damage the internal high-precision electronic components, causing fire or electrical shock.

When plugging or unplugging the USB connector, always hold the connector body. Applying an excessive load to the cable or connector may break the wire or damage this unit.

Handle this unit very carefully.
This unit is assembled with high-precision electronic components. Those components may be damaged if subjected to excessive vibration or shock, causing fire or electrical shock.

About care.
Wipe off dirt on instrument with a soft dry cloth.

1-3. Disposal

When disposing of this unit, take appropriate measures in compliance with applicable regulations regarding waste disposal, and correctly dispose of it yourself, or entrust proper disposal to a licensed industrial waste disposal company. In any case, be sure to comply with the regulations in your country or state to ensure correct disposal.

1-4. Registered trademarks

- Windows 10 Pro are the registered trademarks of Microsoft Corporation in the United States and/or other countries.
- Other product or software names mentioned in this manual may be the trademarks or registered trademarks of respective product manufacturers and a detailed listing is omitted for this manual.

1-5. FCC statement

This device complies with part 15 of the FCC Rules. Operation is subject to the following two conditions: (1) This device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.

Note: This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.

2. Overview


The H17362 series is a photomultiplier tube module with a USB interface. Connecting this unit to a USB port of your PC allows software control to display and record the output data.

The H17362 supports USB hot-swap and so operates without having to restart the PC. The H17362 comes with device drivers and sample software for operation on Windows 10 Pro and so greatly reduces the time required to make setups.

This instruction manual applies to the H17362 series.

3. Part names

The photos below show the external view of the H17362 series.

(1) Light input window

A light-sensitive photocathode is formed on the inner side of this light input window. Do not allow the photocathode to be exposed to excessive light, since excessive light will reduce the sensitivity and service life of this unit.

(2) USB connector

Do not apply excessive shock to the receptacle of this unit when plugging or unplugging the USB connector. Excessive shock may deform or damage the USB connector and the receptacle of this unit. Be careful not to apply excessive force when plugging the USB connector into this unit.

(3) Signal output (SMB receptacle)

Do not apply excessive shock to the receptacle of this unit when plugging or unplugging the SMB receptacle. Excessive shock may deform or damage the SMB connector and the receptacle of this unit. Be careful not to apply excessive force when plugging the SMB connector into this unit.

(4) C mount adaptor (Detachable)

(5) Rod mounting plate (Detachable)

4. Package contents

After unpacking, first check that all items are included. If any part is missing or externally damaged, please notify our sales office.

THz PMT module H17362
C mount adaptor
Rod mounting plate

5. Setup

The procedures explained in this section apply only to cases where your PC runs on Windows 10 Pro.

5-1. Installing the device driver

～The example for the OS Windows 10-Pro～

When using the H17362 first time, install the device driver as described below.

- (1) Turn on the PC power switch.
- (2) Download the device H17362_driver file from following dedicated download page.

URL : https://www.hamamatsu.com/sp/etd/software_etd/download/thz-pmt_module/h17362_en.html

- (3) Copy driver folder and paste it on PC desktop which is stored in (2) downloaded file.
- (4) Connect H17362 to PC with the USB cable of accessory.

Do not connect the USB cable to other than the USB ports.
Connecting the USB cable to other than the USB ports may cause electrical shock or fire, or may damage this unit and PC.

- (5) Right-click the "Start" button in the lower left of the desktop with the mouse and click "Device Manager (M)" from the menu to display the Device Manager.
- (6) Right-click on "HPKUSB_device" in "Other devices" and click "Update Driver(P)" from the menu.
- (7) Click "Browse my computer for H17362_driver software".
- (8) Click the "Browse" button, select the "driver" folder you copied to your desktop, and click "Next".
Installation of the driver software will be started.
- (9) When "The driver has been updated successfully" is displayed, click "Close".

This completes the device driver installation.

5-2. Setting up the sample software

The sample software is stored in the file which downloaded in section 5-1.

Software >> Sample software >> H17362 App.exe

Copy the above executable file and paste it on PC desktop.

See Section 6 “How to use the sample software” for details about the sample software.


The copyright on the sample software is the property of Hamamatsu Photonics. Duplicating and distributing the sample software for any purpose other than controlling the H17362 is not permitted. Adaptation or reverse engineering of this software such as decompilation and disassembly are prohibited.

6. How to use the sample software

6-1. Sample software screen

When you start the sample software that you set up in 5-2 “Setting up the sample software,” the following screen will be appeared.

B1.Read devices button

Reflects the serial number of PMT module to CB1, which is connected to PC.

B2.Set PMT Voltage button

Supply voltage to the PMT module which is indicated in CB1. The voltage is indicated in E1.

B3.PMT OFF button

Stop voltage supply to the PMT module which is indicated in CB1.

CB1.Connected PMT module combo box

Select the serial number of PMT module that to be operated.

The serial number is indicated in CB1 once B1 is clicked.

All devices button designates all connected PMT modules as setting subject.

E1. PMT Voltage[V] edit box

Designates the supply voltage value for PMT module.


D1. Voltage status display

Displays the current PMT voltage supply status ON or OFF.

6-2. Operating procedures

1. Make the necessary setup in advance (See Section 5 “Setup.”)
2. Turn on the PC power switch.
3. Connect the USB cable to the H17362.
Connect the other end of the USB cable to a USB port of the PC.
4. Start “H17362 Sample Software.exe”.
Install H17362 into the experimental setup.

Avoid excessive light because it causes noise and damage.

5. Press “B1Read Devices” button.
6. Select the serial number of H17362 that to be operated in CB1 “Select Port Number combo box”.
7. Enter the HV setting voltage value in the “E1.HV edit box”.
8. Press “B2.HV ON button” to start supply high voltage to the photodetector (PMT).
9. Press “B3.PMT OFF” button to stop the operation of H17362.

Note1. In order to obtain stable operation, it is recommended to perform warm-up operation for several tens of minutes while supplying voltage and injecting THz waves.

Note2. H17362 has polarization characteristics. Please refer 7-4. Photodetector characteristics.

If communication or operation fails during measurement, then restore it by using the following procedure:

1. Open the Task Manager by pressing CTL+ALT+DEL and exit the sample software.
 2. Unplug the USB cable from the H17362.
 3. Plug the USB cable into the H17362.
 4. After the OS has recognized the H17362, restart the sample software.
-

Do not start other application software.
The H17362 is communicating with the PC via USB. If application software that imposes a load on the OS is started, such as those for processing complex calculations or displaying graphs, the desired communication speed may not be obtained, causing the sample software to freeze.

7. Specifications

7-1. General ratings

Parameter	Value	unit
Power input	USB bus power(+4.75 to +5.25)	Vdc
Current consumption	0.1	A
Dimensions	42(W) × 70(H) × 65(D)	mm
Weight	266	g
Location for use	Indoor use	-
Altitude for use	2000 m or less	-
Operating temperature ¹⁾	+5 to +50	°C
Storage temperature ¹⁾	-20 to +50	°C
Storage humidity ¹⁾	85 or less	% RH

1) With no condensation environment

7-2. Specifications

Parameter			Value	unit
Max. supply voltage to PMT			-1800	V
Recommended supply voltage to PMT adjustment range			-1000 to -1800	V
Effective area			□6	mm
Recommended spectral response			0.5 to 2.0	THz
Anode	Threshold E-field ¹⁾	Typ.	5	kV/cm
		Max.	10	
	Polarization dependence		Yes	
Input current ^{1) 3)}			72.5	mA
Settling time ²⁾		Max.	10	s
Support OS			Windows10 Windows11	-

1) Input voltage = +5 V, supply voltage to PMT = -1500 V

2) the time required for the output to reach a stable level following a change in the supply voltage from -1500 V to -1000 V

3) Dark condition

7-3 Other Specifications

Overvoltage category : I 、 Pollution degree : 2、 IP code : IP20

Applicable standards

Safety

IEC 61010-1

EMC

FCC 47CFR Part 15 Subpart B ClassB

IEC 61326-1:


Emission limites:CISPR 11 Group 1, ClassB

Immunity requirements: Table 1



Criterion A	
USB communication	Not occur communication error / stop.
measurement count number	The fluctuation of output must be within -50% to +100 %.

Criterion B	
USB communication	self-recover to normal condition from temporal communication stop USB communication error may occur during test. Once the interference is removed, it must self-recover to the normal condition.
measurement count number	The fluctuation of measurement may temporarily exceed Criterion A during the test. Once the interference is removed, the fluctuation of output must be within -50% to +100%.




Criterion C	
USB communication	USB communication error may occur during test. Once the interference is removed, it must return to normal condition by operator's manual operation.
measurement count number	The fluctuation of measurement may over Criterion A temporarily during the test. Once the interference is removed, it must return to the normal condition by operator's manual operation.

	This mark shown on this unit indicates the power supply is direct current.
---	--

7-4. Photodetector characteristics

H17362 Polarization characteristics


The graph describes relative PMT output while rotating THz wave.

H17362 is designed to output signal for polarized THz waves which parallel to the notch direction in the detection surface.

7-5. Dimensional outlines


(Unit: mm)


CROSS-SECTION VIEW A-A

* Mounting hole for optical rod

7-6. Block diagram

Installation example

In case the optical rod is connected to mounting hole A, the detection part is directly above to the optical rod.

In case the optical rod is connected to mounting hole B, the center position of module body is directly above to the optical rod.

The rod mounting plate is removable. The plate can be attached any plane of H17362 depending on the polarization of interested THz wave.

7-7. DLL file instructions


Outline

The H17362 is a THz PMT module with a built-in USB interface. It is connected to the USB port of your personal computer and controlled by software. By using the API functions described in this instruction manual, you can freely control the PMT applied voltage inside the H17362.

API functions for control and configuration

Software Structure

The figure below shows the configuration of the software that controls this product and acquires count data.

Normal Win32 application software calls H17362_API.dll from within the application.

H17362_API.dll contains API functions starting with “Mid”. This API can be used to control the product. The API is written in C#.

Using Directive

To use each function in C#, declare it as follows.

```
using Middle_SERIAL_NUM; //namespace for H17362 serial number etc. control
using Middle_PMT_VOLTAGE; //namespace for PMT voltage control
using Config; //namespace of common configuration file
```

API Function Overview

The following is a brief description of the API functions of this product.

- | | |
|-----------------------------|---|
| (1)Mid_SN_Get_Connect_Num() | Obtains the number of devices connected to the control PC. |
| (2)Mid_SN_Get() | Obtains the PMT serial number of the H17362 connected to the control PC. |
| (3)Mid_Vcont_Write() | Sets the voltage supply to the light detection element inside the H17362. |

API Function Return Values

The H17362API takes the following enumerations as return values.

Values other than FT_OK are execution errors.

```
public enum FT_STATUS
{
    FT_OK,
    FT_INVALID_HANDLE,
    FT_DEVICE_NOT_FOUND,
    FT_DEVICE_NOT_OPENED,
    FT_IO_ERROR,
    FT_INSUFFICIENT_RESOURCES,
    FT_INVALID_PARAMETER,
    FT_INVALID_BAUD_RATE,
    FT_DEVICE_NOT_OPENED_FOR_ERASE,
    FT_DEVICE_NOT_OPENED_FOR_WRITE,
    FT_FAILED_TO_WRITE_DEVICE,
    FT_EEPROM_READ_FAILED,
    FT_EEPROM_WRITE_FAILED,
    FT_EEPROM_ERASE_FAILED,
    FT_EEPROM_NOT_PRESENT,
    FT_EEPROM_NOT_PROGRAMMED,
    FT_INVALID_ARGS,
    FT_NOT_SUPPORTED,
    FT_OTHER_ERROR,
    FT_TIME_OUT,
    FT_SERIAL_ERROR
};
```

API Function Specifications

(1) Mid_SN_Get_Connect_Num()

<Declaration>

```
FT_STATUS Mid_SN_Get_Connect_Num(ref UInt32 NumDevicesBuf);
```

<Argument>

NumDevicesBuf: Number of device connections buffer

<Return value>

If the function succeeds, the return value is FT_OK.

If the function fails, the return value is anything but FT_OK.

<Explanation>

Check the number of H17362s connected to the control PC.

The number of connections to the control PC is stored in NumDevicesBuf specified in the argument. When controlling multiple H17362s, the NumDevicesBuf specified in the argument is used in Mid_SN_Get() described below.

(2) Mid_SN_Get()

<Declaration>

```
FT_STATUS Mid_SN_Get(UINT32 NumDevices, ref string[] SerialNumsBuf);
```

<Argument>

NumDevicesBuf: Buffer for number of device connections
(1 if H17362 is connected to control PC by itself)
SerialNumsBuf: Serial number storage buffer

<Return value>

If the function succeeds, the return value is FT_OK.
If the function fails, the return value is anything but FT_OK.

<Explanation>

Calls the PMT serial number of the H17362 connected to the control PC.
The PMT serial number of the H17362 connected to the control PC is stored in the SerialNumsBuf specified in the argument.
The number of elements of SerialNumsBuf is required as much as specified in NumDevicesBuf.
When multiple H17362s are connected to multiple PCs, which serial number goes into which element of SerialNumsBuf is determined at random.

(3) Mid_Vcont_Writ()

<Declaration>

```
FT_STATUS Mid_Vcont_Write(UINT32 index, ref string SN, ref UInt16 Vcont);
```

<Argument>

index: ID number of the H17362 to be controlled
(0 if the H17362 is connected to the control PC by itself)
SN: PMT serial number of the H17362 to be controlled
Vcont: PMT control voltage [V] of H17362 (0 to 1800[V])

<Return value>

If the function succeeds, the return value is FT_OK.
If the function fails, the return value is anything but FT_OK.

<Explanation>

Sets the PMT applied voltage inside the H17362 specified by SN, where Index is the device ID randomly assigned when the H17362 is connected to the PC, and matches the SerialNumsBuf element number when Mid_SN_Get() is executed.

8. Warranty

Hamamatsu photomultiplier tube modules and related products are warranted to the original purchaser for a period of one year after delivery. The warranty is limited to repair or replacement of defective products due to defects in workmanship or materials used in their manufacture. Even if within the warranty period, the warranty shall not apply to failures due to misuse, mishandling, modification by the customer, or accidents such as natural or manmade disasters. The customer should inspect and test the products and accessories as soon as they are delivered.

If the equipment is used in a manner not specified by this instruction manual, the protection provided by the equipment may be impaired. We are not liable for damage or losses resulting from failure to comply with precautions or instructions listed in this manual.

9. Repair

If a failure is found in the product or you suspect possible trouble, contact us with the specific symptom and detailed description of the trouble, as well as the production serial number. We will make every effort to repair the returned product in as short a time as possible. In the following cases, however, the repair might be expensive or repair might be refused.

- The product was purchased a long time ago.
- The maintenance parts used in the product are not in current production.
- The product was modified or altered.
- The product was severely damaged.
- Other unavoidable factors.

Product specifications, accessories and contents in this instruction manual are subject to change without prior notice.

HAMAMATSU PHOTONICS K.K. www.hamamatsu.com

Electron Tube Division

314-5, Shimokanzo, Iwata City, Shizuoka Pref., 438-0193, Japan, Telephone: (81)539/62-5248, Fax: (81)539/62-2205

U.S.A.: HAMAMATSU CORPORATION: 360 Foothill Road, Bridgewater, NJ 08807, U.S.A., Telephone: (1)908-231-0960, Fax: (1)908-231-1218

Germany: HAMAMATSU PHOTONICS DEUTSCHLAND GMBH: Arzbergerstr. 10, 82211 Herrsching am Ammersee, Germany, Telephone: (49)8152-375-0, Fax: (49)8152-265-8 E-mail: info@hamamatsu.de

France: HAMAMATSU PHOTONICS FRANCE S.A.R.L.: 19 Rue du Saule Trapu, Parc du Moulin de Massy, 91882 Massy Cedex, France, Telephone: (33)1 69 53 71 00, Fax: (33)1 69 53 71 10 E-mail: infos@hamamatsu.fr

United Kingdom: HAMAMATSU PHOTONICS UK LIMITED: 2 Howard Court, 10 Tewin Road, Welwyn Garden City, Hertfordshire, AL7 1BW, UK, Telephone: (44)1707-294888, Fax: (44)1707-325777 E-mail: info@hamamatsu.co.uk

North Europe: HAMAMATSU PHOTONICS NORDEN AB: Torshamnsgatan 35, 16440 Kista, Sweden, Telephone: (46)8-509-031-00, Fax: (46)8-509-031-01 E-mail: info@hamamatsu.se

Italy: HAMAMATSU PHOTONICS ITALIA S.R.L.: Strada della Moia, 1 int. 6 20044 Arese (Milano), Italy, Telephone: (39)02-93 58 17 33, Fax: (39)02-93 58 17 41 E-mail: info@hamamatsu.it

China: HAMAMATSU PHOTONICS (CHINA) CO., LTD.: 1201, Tower B, Jilaming Center, 27 Dongsanhuan Beilu, Chaoyang District, 100020 Beijing, P.R. China, Telephone: (86)10-6586-6006, Fax: (86)10-6586-2866 E-mail: hpc@hamamatsu.com.cn

Taiwan: HAMAMATSU PHOTONICS TAIWAN CO., LTD.: 13F-1, No.101, Section 2, Gongdao 5th Road, East Dist., Hsinchu City, 300046, Taiwan(R.O.C.) Telephone: (886)3-659-0080, Fax: (886)3-659-0081 E-mail: info@hamamatsu.com.tw