

iPSC-derived Neurons, how are they being used on **FDSS** recently



Stem Cell Research 17 (2016) 306–317

Contents lists available at ScienceDirect

Stem Cell Research

journal homepage: www.elsevier.com/locate/scr

Induced pluripotent stem cell - derived neurons for the study of spinocerebellar ataxia type 3

Susanne K. Hansen^{a,b,*}, Tina C. Stummann^b, Helena Borland^b, Lis F. Hasholt^c, Zeynep Tümer^d, Jørgen E. Nielsen^{c,e}, Mikkel A. Rasmussen^{a,1}, Troels T. Nielsen^e, Justus C.A. Daechsel^b, Karina Fog^b, Poul Hyttel^a

^a Department of Veterinary Clinical and Animal Sciences, University of Copenhagen, Groennegårdsvej 7, 1870 Frø C, Denmark
^b H. Lundbeck A/S, Ørttilavej 9, Valby 2500, Denmark
^c Institute of Cellular and Molecular Medicine, University of Copenhagen, Blegdamsvej 3B, 2200 N, Denmark
^d Applied Human Molecular Genetics, Kennedy Center, Department of Clinical Genetics, Copenhagen University Hospital, Rigshospitalet, Gl. Landevej 7, Glostrup 2600, Denmark
^e Neurogenetics Clinic & Research Laboratory, Danish Dementia Research Centre, Rigshospitalet, University of Copenhagen, Blegdamsvej 9, 2100 Copenhagen, Denmark

Original Research

Functional Characterization of Acetylcholine Receptors Expressed in Human Neurons Differentiated from Hippocampal Neural Stem/Progenitor Cells

Kazuyuki Fukushima^{1,2}, Kazuto Yamazaki¹, Norimasa Miyamoto^{2,3}, and Kohei Sawada^{2,3}

Journal of Biomolecular Screening 1–10
 © 2016 Society for Laboratory Automation and Screening
 DOI: 10.1177/1087057116651166
jbx.sagepub.com
 SAGE

Cellular Dynamics International
 FUJIFILM COMPANY

Spontaneous Ca²⁺ Oscillations in Neuronal Cells

iCell® DopaNeurons – Baseline

Consistent assay signal from well iCell® DopaNeurons are cultured in BrainPhys® neuronal medium (STEM until DIV 14 and then assayed in Mg²⁺- to develop spontaneous Ca²⁺ oscillations

5 μM DNQX

20 μM D-AP5

Calcium bursts are sensitive to receptor blockade. The assay signal is disrupted with antagonists of AMPA receptors, suggesting that these “bursts” are synaptically driven events that are being blocked. Development of this functional endpoint complements other assay technologies (e.g., MEA) used to assess synchronous neuronal activity in human iPSC-derived neuronal culture.

www.nature.com/scientificreports

SCIENTIFIC REPORTS

OPEN **Early pathogenesis of Duchenne muscular dystrophy modelled in patient-derived human induced pluripotent stem cells**

Received: 12 October 2014
 Accepted: 11 May 2015
 Published: 20 August 2015

Emi Shoji^{1,2}, Hidetoshi Sakurai³, Tokiko Nishino⁴, Tatsutoshi Nakahata⁵, Toshio Heike³, Tomonari Awaya³, Nobuharu Fujii⁶, Yasuko Manabe⁵, Masafumi Matsuo⁶ & Atsuko Sehara-Fujisawa²

www.nature.com/scientificreports

SCIENTIFIC REPORTS

OPEN **Calcium dysregulation contributes to neurodegeneration in FTLD patient iPSC-derived neurons**

Received: 21 March 2016
 Accepted: 20 September 2016
 Published: 10 October 2016

Keiko Imamura¹, Naruhiko Sahara², Nicholas M. Kanaan³, Kayoko Tsukita¹, Takayuki Kondo¹, Yumiko Kutoku⁴, Yutaka Ohsawa⁴, Yoshihide Sunada⁶, Koichi Kawakami⁵, Akitsu Hotta¹, Satoshi Yawata⁶, Dai Watanabe⁶, Masato Hasegawa⁷, John Q. Trojanowski⁸, Virginia M.-Y. Lee⁸, Tetsuya Suhara², Makoto Higuchi² & Haruhisa Inoue¹