


Image acquisition software

NDP.scan

Features two “Batch” modes for dealing with large numbers of slides, and a “Single” mode for dealing with individual slides.

The NanoZoomer series comes with its own image acquisition software, NDP.scan, which provides an easy way to convert glass slides into digital ones. You can choose from three different scan modes to fit your applications.



Batch modes

Fully Automatic

Automatically loads and scans each slide.

NDP.scan performs all necessary operations automatically. This includes loading slides, recognizing the sample location, and setting the scan areas, focus points, and other options. It can scan up to the maximum slide capacity (XR: 320 slides, S210: 210 slides, S60: 60 slides) automatically.



Semi Automatic

Automatically scans all slides. You can set specific scan options manually.

NDP.scan will scan the slides automatically after you set the scan options for each slide.

- Configurable Options -
- Scan area
- Scan resolution (20× or 40× mode)
- Number of layers (Z-stack)
- Focus point locations



Single mode

Single Slide Mode

You set scan options for each individual slide.

You can set scan options for each individual slide. The configurable options are the same as those available for Semi Automatic. Also, in Single Slide Mode, you can verify the focal points and make precise adjustments manually to each focus point.





Image acquisition software NDP.scan

Convenient and user-friendly scan features

Profile functions (Saves scanning parameters for each user or sample)




Barcode reading (Automatically reads a slide's barcode information for inclusion in the digital slide)

NDP.scan can automatically read a 1D* or 2D** barcode and include this information in the digital slide. This function makes it easy to share slide data with other databases.


*1D barcode reading is standard; **2D barcode reading is optional.

*1D barcodes: code39, code128
 **2D barcodes: Data Matrix, QRcode, MicroQR
 Please consult with your local sales representative.
 if you would like to inquire about other kinds of barcodes.



Ability to handle samples with uneven surfaces (Adjusts focus as though it were tracing the sample surface)

NDP.scan can adjust the focus as though it were tracing the sample surface. The NanoZoomer series uses a single focal plane for even surfaces, but if the surface of a sample is uneven, the NanoZoomer series will break the sample area up into smaller regions and assign an appropriate focal plane to each area.



* In EU, five types of NanoZoomer (NanoZoomer-XR, NanoZoomer-SQ, NanoZoomer S210, NanoZoomer S60, NanoZoomer S360), NDP.view2 (U12388-21), NDP.view2 Plus (U12388-22) and NDP.serve3 software are CE marked under EU's In Vitro Diagnostics Directive (IVDD) for in vitro diagnostic use.
 In China, three types of NanoZoomer (NanoZoomer 2.0-HT, NanoZoomer 2.0-RS, NanoZoomer-XR) are registered for in vitro diagnostic use.
 In the US, Japan and other countries, NanoZoomer is for research use only and is not permitted to use for clinical diagnostic purposes.

★ NanoZoomer or NDP is a registered trademark of Hamamatsu Photonics K.K. (EU, Japan, U.S.A)

★ Product and software package names noted in this documentation are trademarks or registered trademarks of their respective manufacturers.

★ Subject to local technical requirements and regulations. Availability of products included in this promotional material may vary. Please consult with your local sales representative.

• Information furnished by HAMAMATSU is believed to be reliable. However, no responsibility is assumed for possible inaccuracies or omissions.

• Specifications and external appearance are subject to change without notice.

© 2018 Hamamatsu Photonics K.K.

HAMAMATSU PHOTONICS K.K. www.hamamatsu.com

HAMAMATSU PHOTONICS K.K., Systems Division

812 Joko-cho, Higashi-ku, Hamamatsu City, 431-3196, Japan, Telephone: (81)53-431-0124, Fax: (81)53-433-8031, E-mail: export@sys.hpk.co.jp

U.S.A.: Hamamatsu Corporation, 360 Foothill Road, Bridgewater, NJ 08807, U.S.A., Telephone: (1)908-231-0960, Fax: (1)908-231-1218 E-mail: usa@hamamatsu.com

Germany: Hamamatsu Photonics Deutschland GmbH.: Arzbergerstr. 10, D-82211 Herrsching am Ammersee, Germany, Telephone: (49)8152-375-0, Fax: (49)8152-265-8 E-mail: info@hamamatsu.de

France: Hamamatsu Photonics France S.A.R.L.: 19, Rue du Saule Trapu, Parc du Moulin de Massy, 91882 Massy Cedex, France, Telephone: (33)1 69 53 71 00, Fax: (33)1 69 53 71 10 E-mail: infos@hamamatsu.fr

United Kingdom: Hamamatsu Photonics UK Limited: 2 Howard Court, 10 Tewin Road, Welwyn Garden City, Hertfordshire AL7 1BW, UK, Telephone: (44)1707-294888, Fax: (44)1707-325777 E-mail: info@hamamatsu.co.uk

North Europe: Hamamatsu Photonics Norden AB: Torshamnsgatan 35 16440 Kista, Sweden, Telephone: (46)8-509 031 00, Fax: (46)8-509 031 01 E-mail: info@hamamatsu.se

Italy: Hamamatsu Photonics Italia S.r.l.: Strada della Moia, 1 int. 6, 20020 Arese (Milano), Italy, Telephone: (39)02-93 58 17 33, Fax: (39)02-93 58 17 41 E-mail: info@hamamatsu.it

China: Hamamatsu Photonics (China) Co., Ltd.: 1201 Tower B, Jiaming Center, 27 Dongsanhuang Beltu, Chaoyang District, 100020 Beijing, China, Telephone: (86)10-6586-6006, Fax: (86)10-6586-2866 E-mail: hpc@hamamatsu.com.cn

Taiwan: Hamamatsu Photonics Taiwan Co., Ltd.: 8F-3, No.158, Section2, Gongdao 5th Road, East District, Hsinchu, 300, Taiwan R.O.C. Telephone: (886)03-659-0080, Fax: (886)03-659-0081 E-mail: info@hamamatsu.com.tw

Cat.No.SBIS0064E04
 APR/2018 HPK
 Created in Japan