

## $\beta$ -arrestin recruitment for GPCR Screening

### Introduction

Over the years, it has been demonstrated that signaling through GPCR is not restricted to the type of G protein they are linked to. Downstream events can also be triggered by  $\beta$ -arrestin, a protein involved in signal regulation and receptor internalization. In order to detect  $\beta$ -arrestin recruitment by GPCR activation, DiscoverRx has developed the PathHunter luminescent  $\beta$ -arrestin assay. The latest version of this assay is the PathHunter eXpress  $\beta$ -arrestin protocol which was tested on Hamamatsu FDSS instruments (7000 and 7000EX) and the new FDSS/ $\mu$ CELL. The FDSS/ $\mu$ CELL (a compact, cost effective and simple system for MTS and assay development) and FDSS7000EX are both equipped with Hamamatsu's new EM-CCD camera (Electron Multiplied) for fluorescence and luminescence assays and the FDSS7000 is equipped with our trademark photon counting camera. In this context, we compared FDSS7000 and FDSS7000EX in 1536 format using a proven and reliable assay such as PathHunter. The same experiment was conducted in 384 well format using the FDSS/ $\mu$ CELL. This study shows the comparison (assay window, EC<sub>50</sub> and signal to noise ratio) between the two types of detection systems using the same assay.

### Materials and Methods

One day before the experiment cells were plated in 1536 well clear bottom plates (500 cells/well in 2  $\mu$ l of culture medium) or in 384 well plate (5000 cells per well in 50  $\mu$ l). The day of the experiment, 1  $\mu$ l of compound (3x concentrated diluted in HBSS/20 mM Hepes/0.2 %/fatty acid free BSA) or 5  $\mu$ l in 384 well format (5x concentrated diluted in HBSS/20 mM Hepes/0.2 %/fatty acid free BSA) were added in each well and plates were incubated 90 min at 37 °C.  $\beta$ -arrestin recruitment is detected following the addition of 1  $\mu$ l in 1536 (or 12.5  $\mu$ l in 384 well) of DiscoverRx PathHunter Glow reagent with a multidrop and incubation for 60 min at room temperature. The plates were read on FDSS7000, FDSS7000EX (both in 1536 format) and FDSS/ $\mu$ CELL (in 384 format) with exposure times of 1 s, 5 s and 10 s.






Fig. 1: Curve for PathHunter eXpress  $\beta$ -arrestin GPCR targets run on Hamamatsu FDSS7000 (left) and FDSS7000EX (right) with exposure times of 1, and 10 s

FDSS7000EX

FDSS/μCELL



**Fig. 2:** Dose response curve for PathHunter eXpress  $\beta$ -arrestin GPCR targets run on Hamamatsu FDSS/mCELL with exposure times of 1 s, 5 s and 10 s.

## Results

In this study, we compared the performance of FDSS7000 (photon counting camera) with FDSS7000EX in 1536 format and FDSS/μCELL in 384 format (EM-CCD camera) monitoring the luminescence signal following  $\beta$ -arrestin recruitment. To obtain comparable results, maximum sensitivity was used on the FDSS7000EX, whereas sensitivity 7 (maximum 10) was used on the FDSS7000. Assay window is improved when the assay is run on the FDSS7000EX. No significant difference was observed in EC<sub>50</sub> and S/N ratio between the two FDSS7000 systems. The 1536 format is not implemented on the FDSS/μCELL (made for MTS and assay development) so experiments were conducted in 384. Taking into account the difference in assay format, we found no significant difference between results obtained with the FDSS/μCELL and the FDSS7000.

## Summary

The aim of our study was to compare two different detection methods for luminescence: the FDSS7000 using a photon counting camera, the FDSS7000EX and the FDSS/μCELL, both using an EM-CCD camera for fluorescence and luminescence. We observed no significant difference between the two detection methods and even noticed an improved assay window with the EM-CCD camera. The data shows that even though EM-CCD sensitivity is technically lower than photon counting, luminescent assays such as PathHunter can be successfully run on all our platforms.

★ FDSS is registered trademark of Hamamatsu Photonics K.K. (China, France, Germany, Italy, Japan, U.K., U.S.A.)

★ Product and software package names noted in this documentation are trademarks or registered trademarks of their respective manufacturers.

- Subject to local technical requirements and regulations, availability of products included in this promotional material may vary. Please consult your local sales representative.
- Information furnished by HAMAMATSU is believed to be reliable. However, no responsibility is assumed for possible inaccuracies or omissions. Specifications and external appearance are subject to change without notice.

© 2013 Hamamatsu Photonics K.K.

**HAMAMATSU PHOTONICS K.K.** [www.hamamatsu.com](http://www.hamamatsu.com)

**HAMAMATSU PHOTONICS K.K., Systems Division**

**812 Joko-cho, Higashi-ku, Hamamatsu City, 431-3196, Japan, Telephone: (81)53-431-0124, Fax: (81)53-435-1574, E-mail: export@sys.hpk.co.jp**

**U.S.A.:** Hamamatsu Corporation: 360 Foothill Road, Bridgewater, N.J. 08807-0910, U.S.A., Telephone: (1)908-231-0960, Fax: (1)908-231-1218 E-mail: usa@hamamatsu.com

**Germany:** Hamamatsu Photonics Deutschland GmbH: Arzbergerstr. 10, D-82211 Hensching am Ammersee, Germany, Telephone: (49)8152-375-0, Fax: (49)8152-2658 E-mail: info@hamamatsu.de

**France:** Hamamatsu Photonics France S.A.R.L.: 19, Rue du Saule Trapu, Parc du Moulin de Massy, 91882 Massy Cedex, France, Telephone: (33)1 69 53 71 00, Fax: (33)1 69 53 71 10 E-mail: infos@hamamatsu.fr

**United Kingdom:** Hamamatsu Photonics UK Limited: 2 Howard Court, 10 Tewin Road, Welwyn Garden City, Hertfordshire AL7 1BW, United Kingdom, Telephone: (44)1707-294888, Fax: (44)1707-325777 E-mail: info@hamamatsu.co.uk

**North Europe:** Hamamatsu Photonics Norden AB: Torshamnsgatan 35 SE-164 40 Kista, Sweden, Telephone: (46)8-509-031-00, Fax: (46)8-509-031-01 E-mail: info@hamamatsu.se

**Italy:** Hamamatsu Photonics Italia S.r.l.: Strada della Moia, 1 int. 6, 20020 Arese (Milano), Italy, Telephone: (39)02-93581733, Fax: (39)02-93581741 E-mail: info@hamamatsu.it

**China:** Hamamatsu Photonics (China) Co., Ltd.: B1201 Jiaming Center, No.27 Dongshanhuan Bellu, Chaoyang District, Beijing 100020, China, Telephone: (86)10-6586-6006, Fax: (86)10-6586-2866 E-mail: hpc@hamamatsu.com.cn

Cat. No. SBI0082E02

NOV/2013 IP

Created in Japan (PDF)