

New Tools for Drug Discovery: Multiplexing the monitoring of Ca^{2+} and cAMP signaling in primary cells using the FDSS[®]/ μ CELL.

Introduction

In drug discovery research cell-based assays are a widely accepted tool for high throughput screening of potential drug candidates. So far primary cells are not being used very frequently although they originate from the actual human tissue of interest, genuinely express the relevant drug targets, e.g. cell surface receptors at physiological levels, and carry the endogenous cofactors necessary for efficient and specific signal transduction. Thus, primary cells resemble the in-vivo situation and allow for a higher predictability of the situation in humans.

Here we demonstrate the suitability of primary cells for signaling assays in compound screenings using Clonetics[™] and Poietics[™] primary sensors, i.e. primary and stem cells expressing luminescent biosensor proteins specific for Ca^{2+} and cAMP. These assay formats are used on the FDSS/ μ CELL plate reader to generate pharmacologically relevant data for GPCR ligands.



Materials and Methods

Detecting increases in intracellular concentrations of Ca^{2+} and cAMP

The primary sensors HUVEC calcium biosensor, HMVEC-L calcium biosensor used here transiently express the Ca^{2+} biosensor i-Photina[®] (Axxam) (see www.lonza.com/primarysensors). The GloSensor[™] (Promega) was expressed alone in hMSC (hMSC cAMP biosensor) or co-expressed with i-Photina[®]. Transfection was carried out using the Amaxa[™] Nucleofector[™] technology.

The cells, provided frozen in vials, were thawed and seeded on 96-well or 384-well plates and allowed to recover over night in a humidified cell culture incubator. The i-Photina[®] expressing cells are ready-to-use preloaded with the substrate coelenterazine, GloSensor[™] expressing cells need to be loaded with Promega's GloSensor[™] cAMP reagent prior to the assay.

Agonists and antagonists were bought from Sigma-Aldrich or Calbiochem. Measurement of luminescence was carried out with the FDSS/ μ CELL. Dose-dependent responses and EC/IC₅₀ values were calculated using area under the curve (AUC) integration (GraphPad Prism[®]) for i-Photina[®] or using maximum peak values for GloSensor[™].




Fig. 1: Mechanism of the i-Photina reaction

Incubation of cells, expressing the i-Photina[®] apo-photoprotein, with coelenterazine in the presence of oxygen leads to formation of a stable complex, the active photoprotein. Calcium released from intracellular stores upon stimulation of the cells with agonists via G-protein coupled receptors binds to the photoprotein. The excited photoprotein converts coelenterazine to coelenteramide and emits a flash of blue luminescence.




Fig. 2: GloSensor[®] Technology

Firefly luciferase has been fused to the cAMP-binding domain of human Protein Kinase A (red). Upon binding of cAMP the whole protein molecule undergoes a conformational change. This activates the luciferase domain which converts luciferin into oxyluciferin and emits luminescence. See www.promega.com/glosensor for more information on the live-cell biosensor.

FDSS7000EX

FDSS/μCELL

Settings of the FDSS® /μCELL

Sensitivity was always set to 5, dispense height was 2 mm in all cases, dispense speed was 50 μL/sec for i-Photina® in 384-well and GloSensor™ (96-well only), 150 μL/sec for i-Photina® in 96-well. Detection interval was 2 sec for i-Photina®, 10 sec for GloSensor™.

Optimizing the settings included increasing the dispense speed from 50 to 150 μL/sec for i-Photina® in 96-well plate which led to increase in signal intensity. Altering dispense height in the range from 0.8 mm to 4 mm did not change luminescence intensity.

Slightly decreasing the detection area per well to 5×5 pixels for 384-well and 14×14 pixels for 96-well plates was necessary to reduce luminescence crossover to neighbouring wells on white plates.

Results

Specific stimulation and inhibition of GPCR-mediated signaling

The FDSS®/μCELL has been used to record flash luminescence as well as glow luminescence from different primary sensors in 96-well and 384-well assay formats. The FDSS®/μCELL offers the benefit of dispensing agonists into multiple samples in one step and recording the luminescence of all samples simultaneously.

The primary sensors assay system is well suited to generate pharmacologically relevant data for substances that trigger Ca²⁺-dependent and cAMP-dependent signaling. Agonists for different classes of GPCRs (histamine receptors, purinergic receptors, β₂-adrenergic receptors) clearly showed dose-dependent responses with the HUVEC and HMVEC-L calcium biosensor and the hMSC cAMP biosensor yielding EC₅₀ values consistent with published data (Fig.3 & 4).

Z' factor was well above 0.5 for all three systems (data not shown).

The cAMP-mediated response to isoproterenol, an agonist of the β₂-adrenergic receptor, could be specifically inhibited through propranolol, a known antagonist of this receptor, generating an IC₅₀ (Fig.4) in line with published data.




Fig. 3: HUVEC and HMVEC-L calcium biosensors respond to physiological ligands of cell surface GPCRs.

The assay was performed in a white 384-well plate as described in Materials and Methods. EC₅₀ values (HUVEC: histamine 5.6 μM, ATP 34 μM; HMVEC-L ATP 2.3 μM) are consistent with published data.




Fig. 4: Response of hMSC cAMP biosensor to GPCR ligand can be specifically inhibited by an antagonist.

The assay was carried out in white 96-well plate as described in Materials and Methods. Data calculated from maximum peak values. The IC₅₀ value of 16 nM for propranolol and the EC₅₀ of 65 nM for isoproterenol are in line with published data.

Consecutive detection of calcium and cAMP signaling in one sample

With hMSC co-expressing both biosensors it is possible to monitor Ca^{2+} -dependent and cAMP-dependent signaling in one sample. The cells were triggered simultaneously with two ligands specific for the histamine and for the prostaglandin E_1 receptor. The i-Photina® signal appears a few seconds after addition of the agonists and vanishes before the GloSensor™ luminescence rises which takes several minutes. As the two different luminescence signals appear consecutively they can easily be distinguished.




Fig. 5: Multiplexing detection of intracellular Ca^{2+} and cAMP levels in hMSC co-expressing both biosensors

The assay was carried out in a white 96-well plate as described in materials and methods. Two agonists were dispensed at the same time (indicated by red arrow) into each sample. i-Photina® flash luminescence was recorded immediately every 2 seconds for a total of 36 seconds (A), followed by measuring GloSensor™ glow luminescence every 10 seconds for 25 minutes (B).

Conclusion

Primary cells transiently transfected with either the calcium biosensor i-Photina®, the cAMP-sensitive GloSensor™ or co-transfected with both are available as ready-to-use tools for drug discovery research. They are a groundbreaking robust new assay system for detecting intracellular secondary messenger-dependent signaling upon stimulation with physiological and synthetic agonists. This was demonstrated for four different classes of G-protein coupled receptors in dose-dependent manners for two plate formats. Using the FDSS®/μCELL it has been possible to consecutively monitor on one instrument Ca^{2+} -dependent as well as cAMP-dependent signaling in one sample, simultaneously stimulated with the respective agonists. The FDSS®/μCELL reader is very well suited to reproducibly generate reliable pharmacologically relevant data with low variability from the Clonetics™ and Poietics™ primary sensors assay system.

Authors

Bodo Ortmann¹, Kristin Atze¹, Silke Valentin¹, Nadine Breuer¹, Mathias Kuehn¹, Steffi Franke¹, Pete Stecha², Brock Binkowski², Alex Batchelor¹, Leon de Bruin¹, Jean-Marc D'Angelo³, Thomas Niedereichholz³, Anthony Pitt¹
 1 Lonza Cologne GmbH, Nattermannallee 1, 50829 Cologne, Germany
 2 Promega Corporation, Madison, WI 53711 USA
 3 Hamamatsu Photonics Europe GmbH

FDSS 7000EX

FDSS/μCELL

Legal disclaimer

For research use only. Not for use in diagnostic procedures. i-Photina® is a registered trademark of Axxam (San Raffaele, Italy). Unless otherwise noted, all other trademarks herein are marks of the Lonza Group or its affiliates. The Amaxa™ Nucleofector™ technology is covered by patent and/or patent pending rights owned by Lonza Cologne GmbH. The i-Photina® luminescent calcium biosensor is covered by patent and patent pending rights owned by Axxam. Glosensor™ is a trademark of Promega (Madison, Wisconsin, USA).

★ **FDSS is registered trademark of Hamamatsu Photonics K.K. (China, France, Germany, Italy, Japan, U.K., U.S.A.)**

★ **Product and software package names noted in this documentation are trademarks or registered trademarks of their respective manufacturers.**

- Subject to local technical requirements and regulations, availability of products included in this promotional material may vary. Please consult your local sales representative.
- Information furnished by HAMAMATSU is believed to be reliable. However, no responsibility is assumed for possible inaccuracies or omissions. Specifications and external appearance are subject to change without notice.

© 2013 Hamamatsu Photonics K.K.

HAMAMATSU PHOTONICS K.K. www.hamamatsu.com

HAMAMATSU PHOTONICS K.K., Systems Division

812 Joko-cho, Higashi-ku, Hamamatsu City, 431-3196, Japan, Telephone: (81)53-431-0124, Fax: (81)53-435-1574, E-mail: export@sys.hpk.co.jp

U.S.A.: Hamamatsu Corporation: 360 Foothill Road, Bridgewater, N.J. 08807-0910, U.S.A., Telephone: (1)908-231-0960, Fax: (1)908-231-1218 E-mail: usa@hamamatsu.com

Germany: Hamamatsu Photonics Deutschland GmbH: Arzbergerstr. 10, D-82211 Hensching am Ammersee, Germany, Telephone: (49)8152-375-0, Fax: (49)8152-2658 E-mail: info@hamamatsu.de

France: Hamamatsu Photonics France S.A.R.L.: 19, Rue du Saule Trappu, Parc du Moulin de Massy, 91882 Massy Cedex, France, Telephone: (33)1 69 53 71 00, Fax: (33)1 69 53 71 10 E-mail: infos@hamamatsu.fr

United Kingdom: Hamamatsu Photonics UK Limited: 2 Howard Court, 10 Tewin Road, Welwyn Garden City, Hertfordshire AL7 1BW, United Kingdom, Telephone: (44)1707-294888, Fax: (44)1707-325777 E-mail: info@hamamatsu.co.uk

North Europe: Hamamatsu Photonics Norden AB: Torshamnsgatan 35 SE-164 40 Kista, Sweden, Telephone: (46)8-509-031-00, Fax: (46)8-509-031-01 E-mail: info@hamamatsu.se

Italy: Hamamatsu Photonics Italia S.r.l.: Strada della Moia, 1 int. 6, 20020 Arese (Milano), Italy, Telephone: (39)02-93581733, Fax: (39)02-93581741 E-mail: info@hamamatsu.it

China: Hamamatsu Photonics (China) Co., Ltd.: B1201 JIaming Center, No.27 Dongsanhuan Beilu, Chaoyang District, Beijing 100020, China, Telephone: (86)10-6586-6006, Fax: (86)10-6586-2866 E-mail: hpc@hamamatsu.com.cn

Cat. No. SBIS0086E02
NOV/2013 IP
Created in Japan (PDF)