


A new HTS platform (96-well plate) for early muscle drug discovery

Using the FDSS[®]/ μ CELL EFS from Hamamatsu Photonics, a response to electrical field stimulation of MyoScreen[™] human donors was observed by measuring intracellular calcium transients in real-time. A novel stimulator of myofiber regeneration/repair, recently identified by CYTOO in an “in-house” drug screen was further characterized on the FDSS with MyoScreen[™] system. We found that the hit compound significantly increased levels of myoplasmic calcium compared to non-treated controls upon electrical stimulation. Overall, the Hamamatsu FDSS/ μ CELL electrical stimulation and calcium imaging platform combined with the MyoScreen[™] in vitro myotube model provides a relevant screening system for modulators of calcium levels in the context of skeletal muscle exercise and performance.



<https://cytoo.com/>



MyoScreen[™]


Calcium flux of MyoScreen[™] using the FDSS/ μ CELL 96 ch EFS

MyoScreen[™] responsiveness to electrical stimulation



Captured images of MyoScreen[™]


*Microscopy images offered by CYTOO.




Donor 1 (female)

Donor 2 (female)

Calcium flux recording with EFS (low frequency stimulation)




From twitch to tetanus states



Increasing stimuli frequency results in wave summation and finally complete tetanus reflecting in vivo phenomenon.

Screening of calcium flux modulators Nicardipine: DHPR inhibitor




Calcium flux modulators can be easily assessed in a high throughput way.

FDSS/ μ CELL and EFS option 96-channel electrode array for EFS (electric field stimulation)



▲ FDSS/ μ CELL



▲ 96-electrode array (patent pending)

- Fluorescence / Luminescence measurements while electric field stimulation.
- Variable EFS frequency.
- Variable EFS voltage.

- Stimulate all 96 wells simultaneously
- Cylindrical electrodes
- Change stimulation voltages by column (Patent pending)

* The FDSS/ μ CELL EFS system should not be used for optically detecting/monitoring change in transmembrane potential of the cells.

The FDSS/ μ CELL EFS system should not be used on any cell or cells in which the user or anyone else has expressed target ion channels.

★ FDSS is registered trademark of Hamamatsu Photonics K.K. (China, France, Germany, Italy, Japan, U.K., U.S.A.)

★ Hamamatsu Photonics K.K. and its subsidiaries or representatives do not warrant that the use of our products and those goods does not infringe the intellectual property rights of third parties.

★ Product and software package names noted in this documentation are trademarks or registered trademarks of their respective manufacturers.

★ Subject to local technical requirements and regulations, availability of products included in this promotional material may vary. Please consult your local sales representative.

- Information furnished by HAMAMATSU is believed to be reliable. However, no responsibility is assumed for possible inaccuracies or omissions.
- Specifications and external appearance are subject to change without notice.

© 2019 Hamamatsu Photonics K.K.

HAMAMATSU PHOTONICS K.K. www.hamamatsu.com

HAMAMATSU PHOTONICS K.K., Systems Division

812 Joko-cho, Higashi-ku, Hamamatsu City, 431-3196, Japan, Telephone: (81)53-431-0124, Fax: (81)53-433-8031, E-mail: export@sys.hpk.co.jp

U.S.A.: Hamamatsu Corporation: 360 Foothill Road, Bridgewater, NJ 08807, U.S.A., Telephone: (1)908-231-0960, Fax: (1)908-231-1218 E-mail: usa@hamamatsu.com

Germany: Hamamatsu Photonics Deutschland GmbH: Arzbergerstr. 10, D-82211 Hersching am Ammersee, Germany, Telephone: (49)8152-375-0, Fax: (49)8152-265-8 E-mail: info@hamamatsu.de

France: Hamamatsu Photonics France S.A.R.L.: 19, Rue du Saule Trapu, Parc du Moulin de Massy, 91882 Massy Cedex, France, Telephone: (33)1 69 53 71 00, Fax: (33)1 69 53 71 10 E-mail: infos@hamamatsu.fr

United Kingdom: Hamamatsu Photonics UK Limited: 2 Howard Court, 10 Twin Road, Welwyn Garden City, Hertfordshire AL7 1BW, UK, Telephone: (44)1707-294888, Fax: (44)1707-325777 E-mail: info@hamamatsu.co.uk

North Europe: Hamamatsu Photonics Norden AB: Torshamnsgatan 35 16440 Kista, Sweden, Telephone: (46)8-509 031 00, Fax: (46)8-509 031 01 E-mail: info@hamamatsu.se

Italy: Hamamatsu Photonics Italia S.r.l.: Strada della Moia, 1 int. 6, 20020 Arese (Milano), Italy, Telephone: (39)02-93 58 17 33, Fax: (39)02-93 58 17 41 E-mail: info@hamamatsu.it

China: Hamamatsu Photonics (China) Co., Ltd.: 1201 Tower B, Jiamei Center, 27 Dongsanhuan Beilu, Chaoyang District, 100020 Beijing, P.R.China, Telephone: (86)10-6586-6006, Fax: (86)10-6586-2866 E-mail: hpc@hamamatsu.com.cn

Taiwan: Hamamatsu Photonics Taiwan Co., Ltd.: 8F-3, No.158, Section2, Gongdao 5th Road, East District, Hsinchu, 300, Taiwan R.O.C. Telephone: (886)3-659-0080, Fax: (886)3-659-0081 E-mail: info@hamamatsu.com.tw