

## ■安全上の注意

以下の表に記載されたホローカソードランプには、陰極に表中の元素が含まれております。これらの製品を廃棄する場合は、含有情報を添付のうえ、中間処理・最終処理が可能な処理業者へ委託をお願いいたします。

| 型 名                           | ランプ1本あたりの有害物質使用量 単位 : g (以下) |     |     |     |     |     |     |     |    |    |    |
|-------------------------------|------------------------------|-----|-----|-----|-----|-----|-----|-----|----|----|----|
|                               | C d                          | H g | P b | A s | B e | N i | S b | S e | Cr | Co | Mn |
| L233-11NB(Na)                 |                              |     | 12  |     |     |     |     |     |    |    |    |
| L233-14NU(Si)                 |                              |     | 1   |     |     |     |     |     |    |    |    |
| L233-19NB(K)                  |                              |     | 12  |     |     |     |     |     |    |    |    |
| L233-24NB(Cr)                 |                              |     |     |     |     |     |     |     | 1  |    |    |
| L233-25NU(Mn)                 |                              |     |     |     |     |     |     |     |    |    | 2  |
| L233-27NU(Co)                 |                              |     |     |     |     |     |     |     |    | 1  |    |
| L233-28NQ(Ni)                 |                              |     |     |     |     | 2   |     |     |    |    |    |
| L233-30NQ(Zn)                 |                              |     | 1   |     |     |     |     |     |    |    |    |
| L233-32NU(Ge)                 |                              |     |     |     |     | 1   |     |     |    |    |    |
| L233-33NQ(As)                 |                              |     | 2   | 1   |     |     |     |     |    |    |    |
| L233-34NQ(Se)                 |                              |     | 1   |     |     |     |     | 1   |    |    |    |
| L233-37NB(Rb)                 |                              |     | 11  |     |     |     |     |     |    |    |    |
| L233-38NB(Sr)                 |                              |     |     |     |     | 1   |     |     |    |    |    |
| L233-3NB(Li)                  |                              |     |     |     |     | 1   |     |     |    |    |    |
| L233-48NQ(Cd)                 | 4                            |     | 13  |     |     |     |     |     |    |    |    |
| L233-49NB(In)                 |                              |     | 5   |     |     |     |     |     |    |    |    |
| L233-4NQ(Be)                  |                              |     |     |     | 1   |     |     |     |    |    |    |
| L233-50NQ(Sn)                 |                              |     |     |     |     | 1   |     |     |    |    |    |
| L233-51NQ(Sb)                 |                              |     |     |     |     |     | 1   |     |    |    |    |
| L233-52NQ(Te)                 |                              |     |     |     |     | 1   |     |     |    |    |    |
| L233-55NB(Cs)                 |                              |     | 11  |     |     |     |     |     |    |    |    |
| L233-5NQ(B)                   |                              |     |     |     |     | 1   |     |     |    |    |    |
| L233-80NU(Hg)                 |                              | 1   |     |     |     |     |     |     |    |    |    |
| L233-82NQ(Pb)                 |                              |     | 2   |     |     |     |     |     |    |    |    |
|                               |                              |     |     |     |     |     |     |     |    |    |    |
| L2433-11NB(Na)                |                              |     |     |     |     | 1   |     |     |    |    |    |
| L2433-19NB(K)                 |                              |     |     |     |     | 1   |     |     |    |    |    |
| L2433-24NB(Cr)                |                              |     |     |     |     |     |     |     | 1  |    |    |
| L2433-25NU(Mn)                |                              |     |     |     |     |     |     |     |    |    | 2  |
| L2433-27NU(Co)                |                              |     |     |     |     |     |     |     |    | 1  |    |
| L2433-28NQ(Ni)                |                              |     |     |     |     | 2   |     |     |    |    |    |
| L2433-30NQ(Zn)                |                              |     | 1   |     |     |     |     |     |    |    |    |
| L2433-32NU(Ge)                |                              |     |     |     |     | 1   |     |     |    |    |    |
| L2433-33NQ(As)                |                              |     | 2   | 1   |     |     |     |     |    |    |    |
| L2433-34NQ(Se)                |                              |     | 1   |     |     |     |     | 1   |    |    |    |
| L2433-38NB(Sr)                |                              |     |     |     |     | 1   |     |     |    |    |    |
| L2433-3NB(Li)                 |                              |     |     |     |     | 1   |     |     |    |    |    |
| L2433-48NQ(Cd)                | 4                            |     | 13  |     |     |     |     |     |    |    |    |
| L2433-4NQ(Be)                 |                              |     |     |     | 1   |     |     |     |    |    |    |
| L2433-50NQ(Sn)                |                              |     |     |     |     | 1   |     |     |    |    |    |
| L2433-51NQ(Sb)                |                              |     |     |     |     |     | 1   |     |    |    |    |
| L2433-52NQ(Te)                |                              |     |     |     |     | 1   |     |     |    |    |    |
| L2433-56NB(Ba)                |                              |     |     |     |     | 1   |     |     |    |    |    |
| L2433-5NQ(B)                  |                              |     |     |     |     | 1   |     |     |    |    |    |
| L2433-82NQ(Pb)                |                              |     | 2   |     |     |     |     |     |    |    |    |
|                               |                              |     |     |     |     |     |     |     |    |    |    |
| L733-201NB(Na-K)              |                              |     | 12  |     |     |     |     |     |    |    |    |
| L733-203NU(Si-Al)             |                              |     | 1   |     |     |     |     |     |    |    |    |
| L733-204NQ(Fe-Ni)             |                              |     |     |     |     | 1   |     |     |    |    |    |
| L733-402NQ(CU-CD-PB-ZN)       | 1                            |     | 1   |     |     |     |     |     |    |    |    |
| L733-601NQ(Co-Cr-Cu-Fe-Mn-Ni) |                              |     |     |     |     | 1   |     |     | 1  | 1  | 1  |